# User Manual

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# EAR AND FOREHEAD THERMOMETER

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English ..... 06  
Español .....55

### Product Info

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[www.aldi.us](http://www.aldi.us)

**+ VIDEO**

## QR codes take you where you want to go quickly and easily

Whether you require **product information**, **spare parts** or **accessories**, details on **warranties** or **aftersales services**, or if you want to watch a **product demonstration video** or **360° product view**, our QR codes will take you there in no time at all.

### What is a QR code?

A QR code (QR = Quick Response) is a type of matrix that can be read with a smartphone camera and that contains a link to a website or contact details, for example.

**Advantage:** You do not need to manually enter a website address or contact details.

### How it works

To scan the QR code, all you need is a smartphone with QR code reader software and an internet connection. This type of software can be downloaded for free from your smartphone's app store.

### Try it out now

Just scan the QR code with your smartphone and find out more about the ALDI product you have purchased.

### Your ALDI Service Portal

All details mentioned above can also be found in the ALDI Service Portal at [www.aldi.us](http://www.aldi.us).

Depending on your smartphone plan you maybe charged for the connection.

#### Product Info

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


[www.aldi.us](http://www.aldi.us)

**+ VIDEO**

# Contents

<b>Overview</b> .....	<b>4</b>
<b>Use</b> .....	<b>5</b>
<b>Product contents/device parts</b> .....	<b>6</b>
<b>General information</b> .....	<b>7</b>
Reading and storing the user manual .....	7
Explanation of symbols.....	7
<b>Safety</b> .....	<b>9</b>
Proper use .....	9
Safety notes.....	9
<b>Information on body temperature</b> .....	<b>16</b>
<b>Setup</b> .....	<b>18</b>
<b>Basic functions</b> .....	<b>20</b>
Indicators for body temperature .....	20
Backlight/flashlight function .....	21
Energy-saving mode.....	21
°Celsius/°Fahrenheit setting.....	22
Setting the date and time.....	22
Memory mode.....	25
<b>Operation</b> .....	<b>27</b>
Measuring the body temperature in your ear	27
Measuring forehead temperature.....	29
Measuring object temperature .....	31
<b>Cleaning and maintenance</b> .....	<b>33</b>
<b>Calibration</b> .....	<b>35</b>
<b>Troubleshooting</b> .....	<b>36</b>
<b>Electromagnetic immunity</b> .....	<b>38</b>
<b>Changing the battery</b> .....	<b>46</b>
<b>Technical data</b> .....	<b>50</b>
<b>Disposal</b> .....	<b>53</b>
Disposing of the packaging .....	53
Disposing of old appliances .....	53

## Product contents/device parts

- 1 LCD display
- 2 Backlight (of the LCD display)
- 3 SET button
- 4 Probe tip with infrared sensor
- 5 Push-on adapter for measuring forehead temperature
- 6 Button (👁️) forehead measurement (with backlight)
- 7 Button (👂) ear measurement (with backlight)
- 8 Button (💡) for activating the flashlight
- 9 MEM button
- 10 Signal emission opening
- 11 Battery compartment cover
- 12 Light (flashlight)
- 13 Table rest
- 14 PM indicator
- 15 Memory indicator (📄) (in memory mode)
- 16 Battery indicator (🔋)
- 17 Indicator (📶) for measuring the temperature of objects
- 18 Indicator (👂) for ear measurement
- 19 Indicator (😊) for forehead measurement
- 20 Segment indicator for the measured temperature
- 21 Segment indicator for FE (fever) and CL (clean)

The product contents include two type AAA batteries.

# General information

## Reading and storing the user manual

This user manual accompanies this ear and forehead thermometer DSO 364 (short form “thermometer”). It contains important information on setup and handling.

Before using the thermometer, read the user manual carefully. This particularly applies for the safety instructions. Failure to do so may result in personal injury or damage to the thermometer.

Store the user manual for future use. Make sure to include this user manual when passing the thermometer on to third parties.

## Explanation of symbols

The following symbols and signal words are used in this user manual, on the thermometer or on the packaging.

 **WARNING!**

This signal symbol/word designates a hazard with moderate risk, which may result in death or severe injury if not avoided.

**NOTICE!**

This signal word warns against potential damages to property.

## General information

---

This symbol provides you with useful supplementary information on assembly or operation.

Use and store in a dry state

Batch name:

**LOT** V4916DS0364

Serial number: **SN** 00001  
(sequential number)

Production date: 2016-08 (year, month)

Manufacturer: Handelshaus Dittmann GmbH, Kissinger Straße 68, D-97727 Fuchsstadt/Germany.

Type BF (body float) protection against electrical shock. A type BF application equipment with greater protection against electric shock to the body, but not directly to the heart.

---

# Safety

## Proper use

The thermometer is only designed for measuring the temperature of the human body in the ear or on the forehead as well as for measuring the temperature of liquids (e.g. milk or water) or surfaces. The thermometer is exclusively intended for private use and is not suitable for commercial application.

The patient may be the intended operator. Only use the thermometer as described in this user manual. Any other use is considered improper and may result in damage to property.

The manufacturer or vendor cannot be held liable for damages incurred through improper or incorrect use.

## Safety notes

### **Danger for children and persons with impaired physical, sensory or mental abilities!**

- This thermometer is not intended to be used by persons (including children) with impaired physical, sensory or mental abilities or who lack experience and/or knowledge unless they

are supervised by a person responsible for their safety or have been instructed by such a person on how to use the thermometer.

- Store the thermometer so that it is not accessible for children and young persons under 18 years of age.
- Supervise children to prevent them from playing with the thermometer.
- Ensure that the thermometer does not end up in the hands of children. Children could ingest small parts or batteries and choke. Children could injure themselves when using the thermometer.

### **Risk of injury!**

Improper handling of the thermometer may result in injury.

- The thermometer must not be used in the event of a defect. Do not attempt to modify, disassemble or repair the thermometer or to replace parts of it.
- If you experience any irregularities while using the thermometer, immediately stop using it and consult your physician.

- Do not engage in any other activity when performing a measurement.
- If you are not sure how to use the thermometer, you should consult your physician before using it.
- Keep this user manual for the entire service life of the product and include it if you pass the thermometer on to third parties. Make sure that the user manual is also accessible to third parties. The user manual is part of the thermometer.
- Avoid improper use and use that does not correspond to the designated purpose.
- The thermometer has a class 1 LED (flashlight function). Do not look directly into the LED beam for a prolonged period and do not point the LED beam in the eyes of other persons or animals.
- This thermometer cannot serve as a substitute for medical advice or treatment! The measurements taken are only for reference. If you experience problems with your health, it is imperative that you consult your physician.
- If possible, measure body temperature in the same ear as the tempera-

- ture in the right or left ear may differ.
- Keep the infrared sensor clean and dry and avoid damages (see **fig. C**). This is the only way to obtain accurate measurements.
- Do not touch the infrared sensor and do not blow on it. Dirt on the infrared sensor could result in inaccurate measurement results.
- Do not attempt to measure your body temperature in your ear if your ear is inflamed or infected!
- Do not use the thermometer if it is damaged. Using a damaged device could result in injury, pose serious risks and result in inaccurate measurements.
- Wait a while before measuring temperature if the thermometer was stored in an area with a large temperature difference to the place where you intend to measure temperature. Both the thermometer and persons whose body temperature is being measured should have been in the same environment for at least 30 minutes.
- Do not measure the temperature in your ear if you suffer from a disease

- affecting your ears or if there is still water in your ear, e.g. from showering or swimming.
- Ear wax or a greasy forehead could result in inaccurate measurements. You should therefore clean your ear or forehead before taking your temperature. Then wait approx. 5 – 10 minutes before you perform a temperature measurement.
  - Repeat the measurements approx. every 15 minutes if the readings are unusually low or if you feel they are inaccurate. If possible, use a conventional clinical thermometer to take an oral or rectal temperature reading for comparison. Ask your physician if in doubt.
  - If ear drops or other medications have been administered in one ear, you should measure your body temperature in the other ear.
  - When taking the temperature of persons with hearing aids or ear plugs, it is recommended to wait approx. 30 minutes after they have been removed from the ear before taking their temperature.

**NOTICE!**

**Risk of damage!**

Improper handling of the thermometer may result in damage.

- Do not use any accessories from other devices.
- Do not let the thermometer fall, do not shake it and do not expose it to shocks.
- For accurate readings, ensure in particular that the infrared sensor of the thermometer is protected against damage and exposure to wetness and dirt.
- Do not expose the thermometer to direct sunlight and excessively high or low temperatures. Only use and store the thermometer in a dry environment that is free of contamination.
- Do not open the thermometer and do not attempt to repair it yourself in the event of disturbances or damage. This will void your warranty claim. The thermometer may only be repaired by authorized qualified professionals.
- Do not store the thermometer in areas

where it is exposed to extremely high or low ambient temperatures (under  $-4^{\circ}\text{F}$  ( $-20^{\circ}\text{C}$ ) or over  $131^{\circ}\text{F}$  ( $55^{\circ}\text{C}$ )) or in an extremely dry or moist environment as this may affect the accuracy of the measurements.

- Do not use the thermometer in a damp environment. Keep all liquids away from the thermometer during use and do not in any case immerse the thermometer in liquids.
- Do not use the thermometer if easily flammable substances and gases or explosives are nearby.
- During use, the thermometer may interfere with other electric or electrical equipment. Therefore, do not use it if other electric devices are nearby.
- Do not use a mobile phone near the thermometer when measuring temperature. Please note that portable and mobile HF (high-frequency) communication equipment (e.g. mobile phones) could interfere with medical electrical equipment.
- Do not use the thermometer in the mountains at altitudes of over 9843 feet (3000 meters).

- Do not use the thermometer at distances of less than 5 feet (1.5 meters) to a short-wave or micro-wave device or high-frequency surgical device.
- Medical electrical devices are subject to special precautionary measures with respect to the EMC (electromagnetic compatibility). Therefore, please observe the EMC notes concerning the installation and start-up of the thermometer.

## Information on body temperature

Body temperature refers to the internal temperature of the body (core temperature). It fluctuates over the course of the day and is lowest in the morning and highest in the late afternoon.

Normal body temperature is not reflected by a fixed value, but a range of values.

Age	Normal body temperature in °C	Normal body temperature in °F
0 to 2 years	36.4 – 38.0 °C	97.5 – 100.4 °F
3 to 10 years	36.1 – 37.8 °C	97.0 – 100 °F
11 – 65 years	35.9 – 37.6 °C	96.6 – 99.7 °F
over 65 years	35.8 – 37.5 °C	96.4 – 99.5 °F

The range for normal body temperature differs from person to person and is affected by numerous factors: Age, gender, time of day, place and manner in which the temperature is measured, level of activity, mood and effect of medication and other factors.

Most people have a lower body temperature when they sleep; as the day progresses, their body temperature increases.

Infants and small children have a higher body temperature than adults. As people age, their body temperature decreases. With children, fluctuations in body temperature tend to be higher, briefer and more frequent.

Measuring the temperature at different areas of the body will return different results: The body temperature measured rectally is normally approx. 0.5 °F (0.3 °C) higher than in the mouth and the temperature in the mouth is in turn 0.5 °F (0.3 °C) higher than the body temperature measured in the armpit.

To determine body temperature, it is recommended that you measure your body temperature when healthy on three to four consecutive days at different times of the day. This way, you can determine the normal body temperature, take note of it and, if necessary, compare the temperature readings with the normal readings taken for your body.

### **Determining your in-ear body temperature**


Research has found that the ear is ideal for measuring body temperature as the hypothalamus (the area of the brain that regulates body temperature) and the eardrum are supplied with blood by the same blood vessels. As a result, changes in body

temperature can be determined sooner and more accurately in the ear than other parts of the body. It is therefore beneficial if you measure your body temperature in one of your ears as shown in **fig. B** to the extent possible.

## Setup

The two batteries have already been inserted in the thermometer, but are protected by an insulating strip.

1. Before first use, open the battery compartment cover **11** on the back of the thermometer.
2. Pull the insulating strip out and close the battery compartment cover again.

The thermometer will now emit 4 brief signal tones to indicate that it has switched on and all indicators will initially appear on the LCD display **1** for a short period.

The indicator will then switch to standby mode and display the basic values of time, month

and day of the week as well as current room temperature in alternation. Additionally, the display and key lighting of the buttons (👂) ear measurement and (👤) forehead measurement illuminate for a few seconds.

If the thermometer is switched off, you can switch it on by pressing any button. The thermometer will then switch to standby mode. In this mode, the time, date and room temperature will continually alternate on the display.

To correctly position the display in standby mode and configuration mode, make sure that the front side with the buttons is pointing towards the user. In temperature measuring mode and in memory mode, it should be positioned the other way around: In this case, the display is correctly positioned by holding the thermometer so that the back side is pointing towards the user.

If you would like to immediately switch back to standby mode from configuration or memory mode, you must push the buttons (👂) ear measurement and (👤) forehead measurement at the same time. The thermometer will automatically switch to standby mode if there is no button activity for approx. 1 minute after performing the temperature measurement or, in configuration or memory mode, within approx. 30 seconds.

If “Lo” or “Hi” appears on the display after taking a measurement, the temperature measured has exceeded or underrun the permissible temperature range.

### **IMPORTANT INFORMATION on “Clean warning indicator”:**

If, after performing a number of temperature measurements, the green backlight of the LCD display blinks twice while emitting 2 audible tones and, this is not a malfunction, but a warning indicator to “Clean” the device! For this, read the chapter “‘Clean’ warning indicator”.

## Basic functions

### Indicators for body temperature

The thermometer has a backlight, which immediately illuminates in different colors after measuring your body temperature indicating the range of your body temperature.

**Green backlight:** The body temperature measured is below 37.5 °C (99.5 °F).

**Red backlight:** The body temperature measured is within the range of an elevated temperature of 37.5 °C (99.5 °F) to 37.9 °C (100.3 °F).

**Red backlight red and additionally FE (fever indicator) on the display:**

The body temperature measured is within the fever range between 38.0 °C (100.4 °F) and 42.9 °C (109.3 °F).

## Backlight/flashlight function

**Backlight:** Every time the SET button **3** or the MEM button **9** is pushed, the backlight **2** of the LCD display turns on for approx. 10 seconds.

**Flashlight:** Pushing the button **8** (💡) for approx. 2 seconds in standby mode, the flashlight **12** of the thermometer (1 LED next to the probe tip) will illuminate for approx. 10 seconds. At the same time, the blue backlight of the buttons (👤) forehead measurement **6** and (👂) ear measurement **7** illuminate for approx. 15 seconds.

If you briefly push the button (💡) in standby mode, the key lighting of the buttons (👂) ear measurement and (👤) forehead measurement illuminate for approx. 15 seconds.

If the thermometer is switched off and you push the button (💡), the backlight of the LCD display will illuminate for approx. 10 second and the backlights of the buttons (👤) forehead measurement and (👂) ear measurement will illuminate for approx. 15 seconds.

## Energy-saving mode

By default, the thermometer is configured to switch to standby mode approx. 30 – 60 seconds after the last button activity depending on the respective mode (time, date and room temperature are shown in alternation). If you would like the thermometer to turn off completely after approx. 1 minute, follow the steps below:

1. In standby mode, hold the SET button **3** down for approx. 2 seconds until the 24 hr indicator starts to blink on the display.
2. Then push the SET button 6 more times until the LCD indicator briefly appears on the LCD display and until the ON indicator flashes right afterwards.
3. If you now push the MEM button **9**, the OFF indicator will start to blink on the display.
4. Then push the SET button to save this setting. The thermometer will now return to standby mode and will now automatically switch off approx. 1 minute after the last button activity.
5. If you would like to set the thermometer back to continuously standby mode, you can set the LCD/OFF setting to ON according to the description above.

## °Celsius/°Fahrenheit setting

1. Push the MEM button **9** to switch to memory mode.
2. Then hold the SET button **3** down for approx. 2 seconds until the temperature measurement/indicator switches from °C (degrees Celsius) to °F (degrees Fahrenheit/1 signal tone) or vice-versa from °F (degrees Fahrenheit) to °C (degrees Celsius/2 signal tones).

## Setting the date and time

In standby mode and in memory mode, the thermometer displays the room temperature, time and date in alternation. The two blinking dots between the hours and minutes indicate the second intervals.

1. Push the SET button **3** down for approx. 2 seconds to switch to the configuration mode. The indicator for the time mode 12 hour (12 hr) or 24 hour (24 hr) will blink.
2. Push the MEM button **9** to select the desired time mode. Push the SET button to save the selected time mode. If the 12 hour indicator (12 hr) is configured, the symbol PM will appear on the display when displaying the time.
3. The indicator for hours will blink on the display and with the MEM button, you can perform the setting from 1 – 12 or 1 – 24 depending on the previously selected time mode. You can then save the selected value with the SET button.
4. Now, the indicator for minutes will blink on the display and with the MEM button, you can set a value from 00 – 59 and save the selected value with the SET button.
5. Now, the indicator for the calendar year will blink on the display and with the MEM button, you can set a value from 2001 – 2099 and save the selected value with the SET button.
6. Now, the indicator for the month will blink on the display and with the MEM button, you can set a value from 1 – 12 and save the selected value with the SET button.
7. Now, the indicator for the weekday will blink on the display and with the MEM button, you can set a value from 1 – 31 and save the selected value with the SET button.

8. Finally, the LCD indicator appears on the display for a short time, immediately followed by a blinking ON. With the MEM button, you can choose between ON (thermometer stays in continuously standby mode) and OFF (thermometer automatically switches off approx. 1 minute after the last button activity) and save the selected setting with the SET button. The thermometer will then switch to standby mode. If ON was selected, it will remain in standby mode continuously. If OFF was selected, it will always switch off automatically after approx. 1 minute of button inactivity.

If you push the MEM button when setting the values for the date and time, the next higher value will be set in each case. Holding the MEM button down will quickly cycle through the configurable values. Releasing the MEM button will leave the setting at the respective value reached.


If you briefly push both buttons (👂) ear measurement and (👤) forehead measurement at the same time in memory or configuration mode, the thermometer will switch back to standby mode and the time, date and room temperature will be shown in alternation. It will also automatically switch to standby mode after approx. 30 seconds of button inactivity.

## Memory mode

The thermometer saves the last 9 body temperature measurements in your ear or on your forehead. The current measurement will appear last in the memory. If all 9 memory slots are occupied, the oldest measurement will be deleted from the memory.

1. In standby mode, hold the thermometer so that the back of the thermometer points towards you to display the saved measurements and push the MEM button **9**.

The memory indicator (🗄️) **15** (memory access) and the number for the last measurement saved (e.g. 8) will now briefly appear on the display.


Directly afterwards, the temperature measured with the indicator for ear measurement (👂) **18** or forehead measurement (😊) **19** will be displayed,

followed by the time

and then the date in continuous alternation.

2. Each time you push the MEM button again, first the number and then the data saved for the next memory slot will be shown briefly in descending order. Once the oldest memory slot with the number 1 is reached, the indicator will switch back to standby mode after pushing the MEM button again and back to the most recent memory slot (e.g. 8) if the MEM button is pushed again etc.

If no measurements have been saved, (---°C) with the memory indicator () will appear on the display.

The following symbols may appear on the display when accessing the measurements in the memory:

Symbol	Explanation
	Memory indicator for memory mode/ memory access
<b>PM</b>	PM indicator (“post meridiem” = afternoon in 12-hour mode)
	Saved ear temperature measurement

Saved forehead temperature measurement


If you push the SET button down for approx. 2 seconds while the saved measurement data is shown, the respective temperature in the memory will switch/ be converted from °C (degrees Celsius) to °F (degrees Fahrenheit) or vice-versa. The switch from °C to °F is confirmed by 1 short audible tone and the switch from °F to °C is confirmed by 2 short audible tones.

## Operation

### Measuring the body temperature in your ear


1. Before measuring your body temperature in your ear, remove the push-on adapter for measuring forehead temperature **5** from the thermometer.
2. Before performing the measurement, make sure that the infrared sensor **4** is clean (see **fig. C**).
3. Hold the thermometer with your thumb or index finger on the button **9** ear measurement **7**, but do not push the button **9** yet.


4. Pull your ear back slightly to straighten your ear canal.

5. Slowly and carefully insert the probe tip with the infrared sensor **4** in your outer ear canal.

6. Hold the thermometer so that the probe tip is directly aligned with the eardrum inside your ear (see **fig. B**). If possible, always perform the temperature measurement in the same manner on the same ear.
7. Hold the thermometer steady and then push the button (👂) ear measurement **7**.

A signal tone will indicate that the measurement has started. A signal tone will sound twice after approx. 2 seconds indicating that the temperature measurement is complete. The indicator (👂) **18** and the measured ear temperature will appear on the LCD display **1**.

At the same time, the backlight **2** and the illuminated ring will light up indicating the range of the measured temperature; see chapter “Indicators for body temperature”.


- Carefully remove the probe tip of the thermometer from your outer ear canal and read the measured temperature on the LCD display. The indicator (  ) indicates that the temperature in your ear has been measured. It is only possible to perform a new measurement after the backlight has gone out (you must wait approx. 6 seconds).

Wait until the measurement has been saved and no circles are visible on the display.

## Measuring forehead temperature


- Fit the push-on adapter for measuring forehead temperature **5** on the thermometer so that it is secure.
- Before performing the measurement, make sure that the infrared sensor is clean. Hold the thermometer with your thumb on the button (  ) forehead measurement **6**, but do not push the button (  ) yet.


## Operation

---

3. Hold the thermometer in the middle of your forehead so that the push-on adapter for measuring your forehead temperature is gently resting on your forehead.

4. Hold the thermometer steady and then push the button (☛) forehead measurement.

A signal tone will indicate that the measurement has started. A signal tone will sound twice after approx. 2 seconds indicating that the temperature measurement is complete. The indicator (☺) 19 and the measured forehead temperature will appear on the display 1.

At the same time, the backlight 2 will light up indicating the range for the measured temperature; see chapter "Indicators for body temperature".

- Carefully remove the probe tip **4** of the thermometer from your forehead and read the measured temperature on the LCD display.

The symbol (☺) **19** indicates that your forehead temperature has been measured. You can now immediately start a new measurement.


If you want to immediately switch back to standby mode, you can briefly push the MEM button. The thermometer will also automatically switch back to standby mode after approx. 1 minute of button inactivity.

## Measuring object temperature

- Please note: This method of performing measurements is not suitable for measuring body temperature!
- Before measuring the temperature of an object (liquid/surface), remove the push-on adapter for measuring forehead temperature from the thermometer.
- Before performing the measurement, make sure that the infrared sensor is clean.
- Push the buttons (👂) ear measurement **7** and (👁️) forehead measurement **6** at the same time. The indicator (🔒) for measuring the temperature of objects **17** accompanied by the indicator “food” on the display shows that the object temperature measurement has been activated.

5. Hold the probe tip with the infrared sensor **4** so that it can measure the temperature of the liquid or surface within a distance of 0.39" (1 cm).

6. Now push the button **7** ear measurement. A signal tone will indicate that the measurement has started. A signal tone will sound twice after approx. 2 seconds indicating that the temperature measurement is complete. The object temperature measured will appear on the display.

If you would like to immediately switch back to standby mode, you can push the buttons **7** ear measurement and **8** forehead measurement at the same time. This will also be done automatically after approx. 1 minute of button inactivity.

## Cleaning and maintenance

### NOTICE!

#### Risk of short circuit!

Water that has penetrated the housing may cause a short circuit.

- Never submerge the thermometer in water.
- Make sure that no water penetrates the housing.


### NOTICE!

#### Risk of damage!


Improper handling of the thermometer may result in damage.

- Do not use any aggressive cleaners, brushes with metal or nylon bristles, sharp or metallic cleaning utensils such as knives, hard scrapers and the like. They could damage the surfaces.
  - Never attempt to clean the thermometer with a scouring agent, petrol or a solvent.
1. For accurate measurements, it is very important that the infrared sensor **4** of the thermometer is clean and free of scratches.

2. Each time before performing a measurement or if the “Clean” warning indicator appears on the display, check the infrared sensor and, if necessary, carefully clean it with a cotton swab soaked with a small amount of medical alcohol and then immediately with a clean, dry cotton swab.

3. Use a soft, dry cloth to clean the housing.


4. After cleaning, store the thermometer in a clean, dry place at room temperature or set it on the table rest **13**.
5. Never subject the thermometer to extreme temperatures, moisture, direct sunlight or vibrations!

### **“Clean” warning indicator**

The thermometer has a special warning indicator that reminds the user when to clean the infrared sensor on the probe tip of the thermometer.

If you have not used the thermometer for a prolonged period of time or after approx. 10 measurements, the warning notice will appear by the segment indicator for CL **21** on the top left of the LCD display **1** and the green backlight **2** of the thermometer will blink twice and emit a signal tone twice to remind the user that the probe tip with the infrared sensor should be cleaned.

## Calibration

The thermometer was calibrated in the factory following manufacture. If the thermometer is used in accordance with the provisions and notices in the user manual for private and non-commercial use, there is no need to calibrate it at regular intervals. If, however, you doubt the accuracy of the measurements, you should immediately contact your vendor.

Never attempt to disassemble the thermometer or make changes as they could result in risks and in such a case, also void any warranty claim!

# Troubleshooting

Fault	Cause	Solution
The batteries have been inserted, but the LED display <b>1</b> is blank.	There may be debris in the battery compartment.	Remove the debris from the battery compartment.
	The batteries are empty or were inserted with incorrect polarity.	Replace the batteries with full batteries. Make sure the polarity is correct!
	The battery contacts are not aligned with the terminals.	Check the battery contacts.
	There is a fault in the electronics.	Remove the batteries and put them back in after approx. 3 seconds.
 <p>The image shows a dashed rectangular box containing two rows of LED display characters. The top row shows 'Lo' and the bottom row shows 'Hi'. Each row has a small circular icon to its left.</p>	The temperature measured is less than 32.0 °C/89.6 °F or higher than 42.9 °C/109.3 °F.	Observe the notices in the user manual.
 <p>The image shows a dashed rectangular box containing two rows of LED display characters. The top row shows 'ErrP' and the bottom row shows 'ErrH'. Each row has a small circular icon to its left.</p>	The thermometer is defective (hardware error).	Contact the manufacturer/vendor.

Fault	Cause	Solution
	The operating temperature is outside of the valid range.	Observe the operating conditions (see chapter “Technical data”).
The temperature measured is not accurate or you doubt the accuracy of the measurement.	Check whether the infrared sensor  is dirty.	Clean the infrared sensor (see chapter “Cleaning and maintenance”).
	Check whether the thermometer is properly used.	Observe the notices on operating the thermometer and obtaining a correct measurement.
	Check if the device and patient have been in the same room for at least 30 minutes before performing the measurement.	Make sure that the device and patient have been in the same room for at least 30 minutes before performing the measurement.
	The thermometer is being used outdoors.	Only use the thermometer in indoor areas.
	You have held the thermometer in your hands for a prolonged period.	Place the thermometer e.g. on a table so that it can adjust to the room temperature.

# Electromagnetic immunity

Table 1 – Instructions and manufacturer data – electromagnetic emissions – for all EQUIPMENT and SYSTEMS (see 5.2.2.1).

## Instructions and manufacturer data – electromagnetic emissions

The (EQUIPMENT or SYSTEM) is designed for use in the electromagnetic environment described below. The customer or user of the (EQUIPMENT or SYSTEM) should ensure that it is indeed used in such an environment.

Emission test	Compliance	Electromagnetic environment – instructions
HF emissions CISPR 11	Group 1	The (EQUIPMENT or SYSTEM) only uses HF energy for internal operation. Therefore, very limited HF emissions may result and it is highly unlikely that they will disrupt nearby electronic equipment in any way.

## Guidelines and manufacturer's declaration – electromagnetic emissions

The DSO 364 model is intended for operation in an environment as described below. The customer or the user of the DSO 364 model should ensure that it is operated in such an environment.

Measuring emitted interference	Compliance	Electromagnetic environment – guidelines
HF emissions in accordance with CISPR 11	Group 1	
HF emissions in accordance with CISPR 11	Class B	The DSO 364 model is intended for use in all facilities including residential areas and for areas that are directly connected to a public utility network that also supplies buildings used for residential purposes.
Harmonic current emissions in accordance with IEC 61000-3-2	Not applicable	
Emissions of voltage fluctuations/flicker in accordance with IEC 61000-3-3	Not applicable	

### **Guidelines and manufacturer's declaration – electromagnetic immunity**

The DSO 364 model is intended for operation in an electromagnetic environment as described below. The customer or the user of the DSO 364 model should ensure that it is used in such an environment.

## Electromagnetic immunity

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic discharge (ESD) in accordance with IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	Not applicable ± 8 kV air discharge	Floors should be made of wood or concrete or be covered with ceramic tiles. If the floor is covered with synthetic material, the relative humidity must amount to at least 30 %.
Fast transient electrical bursts in accordance with IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines	Not applicable	The quality of the supply voltage should correspond with a typical store or hospital environment.
Surges in accordance with IEC 61000-4-5	± 1 kV symmetrical voltage ± 2 kV asymmetrical voltage	Not applicable	The quality of the supply voltage should correspond with a typical store or hospital environment.

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Voltage drops, short interruptions and fluctuations in the supply voltage in accordance with IEC 61000-4-11	$< 5\% U_T$ ( $> 95\%$ drop of $U_T$ ) for $\frac{1}{2}$ period $40\% U_T$ ( $60\%$ drop of the $U_T$ ) for 5 periods $70\% U_T$ ( $30\%$ drop of $U_T$ ) for 25 periods $< 5\% U_T$ ( $> 95\%$ drop in $U_T$ ) for 5 s	Not applicable	The quality of the supply voltage should correspond with a typical store or hospital environment. If the user of the DSO 364 model also requires that the function continues after interruptions to the power supply, it is recommended that you supply the DSO 364 model from an uninterrupted power supply or battery.
Magnetic field with supply frequency (50/60 Hz) in accordance with IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields associated with the mains frequency should correspond to the typical levels encountered in a store or hospital environment.
<b>NOTICE:</b> $U_T$ is the mains AC voltage prior to application of the test levels.			

### Guidelines and manufacturer's declaration – electromagnetic immunity

The DSO 364 model is intended for operation in an electromagnetic environment as described below. The customer or the user of the model should ensure that it is used in such an environment.

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
			<p>Portable and mobile wireless equipment should not be used at distances to the [device or system] including its lines that are below the recommended safety clearance that is calculated based on the equation that applies for the transmission frequency. <b>Recommended safety clearance:</b></p>

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Conducted HF disturbances in accordance with IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d = 1.2 \sqrt{P}$
Emitted HF disturbances in accordance with IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz
			<p>where P is the nominal output of the transmitter in watts (W) according to the specifications of the transmitter manufacturer and d is the recommended safety clearance in meters (m).</p> <p>According to an on-site survey<sup>a</sup>, the field strength</p>

			<p>of stationary radio transmitters is less than the compliance level<sup>b</sup> for all frequencies. Disturbances are possible within the vicinity of devices with the following graphic symbols.</p>
--	--	--	---

**NOTICE 1** In the case of 80 MHz and 800 MHz, the higher value shall apply.

**NOTICE 2** These guidelines may not apply in all situations. The propagation of electromagnetic waves is affected by absorption and reflection by buildings, objects and people.

**a.** The field strength of stationary transmitters such as base stations for cordless phones and land mobile radio services, amateur stations, AM and FM radio and television stations cannot be precisely predetermined in a theoretical respect. It is recommended that the location be surveyed in order to determine the electromagnetic environment based on stationary HF transmitters. If the field strength at the location of the DSO 364 model exceeds the compliance level specified above, the DSO 364 model must be observed with respect to normal operation at each application location. If the device behaves unusually, further measures may have to be taken, e.g. re-orientation or relocation of the DSO 364 model.

**b.** Not applicable over a frequency range from 150 kHz to 80 MHz.

Recommended safety clearance between portable and mobile HF telecommunications devices and the [DEVICE or SYSTEM].

The DSO 364 model is intended for operation in an electromagnetic environment where HF disturbances are monitored. The customer or user of the DSO 364 model can help to avoid electromagnetic disturbances by observing the minimum clearance between portable and mobile HF telecommunication devices (transmitters) and the DSO 364 model – depending on the output power of the communication device as discussed below.

Nominal output of the W transmitter	Safety clearance depending on the transmission frequency m		
	150 kHz to 80 MHz $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d=2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters whose nominal output is not specified in the above table, the distance can be determined using the equation in the respective column where P represents the nominal output of the transmitter in watts (W) as specified by the transmitter manufacturer.

**NOTICE 1** In order to calculate the recommended safety clearance between transmitters within a frequency range of 80 MHz to 2.5 GHz, an additional factor of 10/3 was used to reduce the probability that a mobile/portable communication device accidentally brought into the patient area causes a disturbance.

**NOTICE 2** These guidelines may not apply in all situations. The propagation of electromagnetic waves is affected by absorption and reflection by buildings, objects and people.

## Changing the battery

### **Danger of explosion!**

Improper handling of the batteries may cause an explosion.

- Only replace the batteries with the same or an equivalent battery type.
- Do not expose the batteries to excessive heat such as from sunshine or fire.
- The enclosed batteries must not be charged or reactivated with other means, they must not be dismantled, disposed of in fire or short circuited.
- Protect batteries against mechanical damage.

 **WARNING!**

**Risk of injury!**

Handling batteries improperly may result in injury.

- Store the batteries so that they are not accessible for children.
- If swallowed, batteries could pose a risk of death. Therefore, store the batteries so that they are not accessible for small children. If a battery is swallowed, you must seek immediate medical attention.
- If battery acid has leaked, avoid contact with your skin, eyes and mucous membranes. In the event of contact with acid, immediately flush the applicable areas with plenty of clean water and promptly consult a physician.

**NOTICE!**

**Risk of damage!**

Improper handling of the batteries may result in damage.

- Make sure the polarity is correct when inserting the batteries.

- Only insert batteries of the same type; do not use different types of batteries or old batteries with new batteries.
- Remove the old batteries from the battery compartment if they are empty or if you do not intend to use the thermometer for a prolonged period of time.

If the battery indicator (  16 appears on the LCD display  or the thermometer no longer functions, the batteries should be replaced immediately.


All settings and entries in the memory will be erased if the batteries are empty or when changing weak batteries! Therefore, after inserting new batteries, the date, time and energy-saving mode must be reset; see chapter “Setting the date and time”.

1. Type AAA alkaline batteries are required for the thermometer. Do not use any rechargeable batteries!
2. To open the battery compartment cover , push against the locking clip on the back of the thermometer and take the cover off.
3. Remove the empty batteries.
4. Then insert two new alkaline batteries (type AAA). When inserting the batteries, make sure that the polarity of the batteries is correct (see marking/stamped information in the battery compartment).

5. Put the battery compartment cover back in and push it against the thermometer housing until it locks back into place.

## Technical data

Model/type:	DSO 364
Dimensions (L × W × H):	Approx. 5.37" x 1.75" x 1.59" (136.5 × 44.5 × 40.5 mm)
Weight:	Approx. 0.12 lb (53.5 g) (without batteries and table rest)
Material:	Plastics, metals
Power supply:	3.0 V DC, 2× AAA batteries (V = volts, DC = direct current)
Measurement method:	Contactless infrared measurement
Measurement time:	Approx. 2 seconds
Measuring distance:	Approx. 0.39" (1 cm) (mea- surement of objects)
Display:	LCD display with green backlight, data shown in °C (degrees Celsius) or °F (degrees Fahrenheit), res- olution: 0.1 °C/°F
Memory function: (memory)	Storage of the last 9 body temperature measurements
Energy-saving function:	The thermometer au- tomatically deactivates after approx. 1 min. if energy-saving mode is activated.
Laser class:	The LED light complies with protection class I (integrat- ed flashlight)

**Body temperature:**

Measurement range:	32.0 °C – 42.9 °C (89.6 °F – 109.3 °F)
Measurement accuracy:	from 35.5 °C – 42.0 °C (95.9 °F – 107.6 °F): +/-0.2 °C (+/-0.4 °F) other ranges: +/-0.3 °C (+/-0.5 °F)

**Object temperature:**

Measurement range:	20.0 °C – 60.0 °C (68.0 °F – 140.0 °F)
Measurement accuracy:	+/-1.0 °C (+/-2.0 °F)


**Room temperature:**

Measurement range:	5.0 °C – 59.9 °C (41.0 °F – 139.8 °F)
Measurement accuracy:	+/-1.0 °C (+/-2.0 °F)

**Operating conditions:**

Atmospheric pressure:	700 hPa – 1060 hPa
Forehead measurement:	15.0 °C – 40.0 °C (59.0 °F – 104.0 °F) at 20 % – 85 % relative humidity
Ear/object measurement:	10.0 °C – 40.0 °C (50.0 °F – 104.0 °F) at 20 % – 85 % relative humidity

### Storage/transport data:

Storage/transport temperature:  $-20\text{ }^{\circ}\text{C}$  –  $55\text{ }^{\circ}\text{C}$   
at 20 % – 90 % relative humidity

Atmospheric pressure:  
700 hPa – 1060 hPa

### Device nameplate:

Ear-and-Forehead Thermometer DSO 364

 3V DC (2xAAA)

Made in China

**LOT** V4916DSO364 item #94078

Manufactured for:

Handelshaus Dittmann GmbH

Kissinger Straße 68

D-97727 Fuchsstadt/

Germany

**3 Year Warranty**

2016-08

Keep dry

Distributed by ALDI Inc.

Batavia, IL 60510

Visit us at [www.Aldi.us](http://www.Aldi.us)

---

# Disposal

## Disposing of the packaging

- Sort the packaging before you dispose of it. Dispose of paperboard and cardboard with the recycled paper service and wrappers with appropriate collection service.

## Disposing of old appliances

- Should the thermometer no longer be capable of being used at some point in time, **dispose of it in accordance with the regulations in force in your city or state.**
- Please ensure your recycling information applies to local regulations and the EPA recommendations ([www.epa.gov](http://www.epa.gov)).

### **Batteries may not be disposed of with household waste!**

As the end user you are required by law to bring all batteries and rechargeable batteries, regardless whether they contain harmful substances\* or not, to a collection point run by the city, county, or state, or to a retailer, so that they can be disposed of in an environmentally friendly manner.

\*labeled with: Cd = cadmium, Hg = mercury,  
Pb = lead

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