User Manual Manual del usuario


BLOOD PRESSURE MONITOR

TENSIÓMETRO

English 03
Español 63

Contents

- | | | | |
|----------|------------------------|----------|-------------|
| 1 | Blood Pressure Monitor | 3 | Storage bag |
| 2 | Battery (1.5 V/AA), 4× | 4 | Cuff |

ACTIV ENERGY

The delivery contents include four 1.5 V AA batteries.

Table of Contents

Contents	2
General information	5
Reading and storing the user manual	5
Proper use	5
Explanation of symbols	6
Safety	7
Explanation of notes.....	7
General safety information	8
Battery notes.....	19
Safe use of the power supply unit.....	22
Product description	24
Description of display	26
First use	27
Checking the blood pressure monitor and package contents.....	27
Operation	28
Inserting or changing batteries	28
Using a power supply unit.....	30
System settings	30
Notes on use	31
Performing a measurement	34
Fitting the cuff	35
Assuming the correct posture	36
Starting a blood pressure measurement.....	37
Deleting or storing measurement results.....	38
Deactivation	38

Evaluating the measurements	39
About blood pressure.....	39
Cardiac arrhythmia.....	40
Average of measurement results.....	41
Checking memory.....	41
Deleting memory.....	42
WHO classifications.....	42
Questions and answers	43
Troubleshooting	45
FCC Compliance Statement	46
Information about electromagnetic immunity	47
Cleaning and maintenance	53
Cleaning.....	54
Cleaning the cuff.....	54
Changing the batteries.....	54
Storage.....	54
Technical data	56
Full display reference	59
Manufacturer information	59
Disposal	60
Disposing of the packaging.....	60
Disposing of the blood pressure monitor.....	60

General information

Reading and storing the user manual

This user manual accompanies this blood pressure monitor. It contains important information about using and handling the product.

Before using the blood pressure monitor, read the user manual carefully. This particularly applies to the safety instructions. Failure to heed this user manual may result in severe injury or damage to the blood pressure monitor.

Store the user manual for future use. Make sure to include this user manual when passing the blood pressure monitor on to third parties. When passing the device on to third parties, also advise them of the risks that can arise in connection with the blood pressure monitor.

Proper use

The blood pressure monitor is only designed for non-invasive measurement of systolic, diastolic blood pressure and pulse rate of adults, using the oscillometric method. The results are shown on the display of the blood pressure monitor.

The blood pressure monitor issues a warning in the event of possible cardiac arrhythmia. The levels measured are classified and rendered in accordance with the guidelines of the WHO (World Health Organization).

The blood pressure measurement taken with this blood pressure monitor corresponds to the one of a trained observer according to the cuff/stethoscope auscultation method within the limits defined by the Recognized Consensus Standard for electronic sphygmomanometers.

The blood pressure monitor is not intended to be used for toddlers and children.

Explanation of symbols

No special training or knowledge is required to use the blood pressure monitor.

The blood pressure monitor is exclusively intended for private use and is not suitable for commercial or clinical use. The measurement position is on an adult's upper arm only.

Only use the blood pressure monitor as described in this user manual. Any other use is deemed improper and may result in damage to property or even personal injury. The blood pressure monitor is not a children's toy.

The manufacturer or vendor accepts no liability for damage caused by improper or incorrect use.

Explanation of symbols

The following symbols and signal words are used in this user manual, on the blood pressure monitor or on the packaging.

This symbol provides you with useful supplementary information about use.

This symbol identifies devices classified as type BF.


This symbol indicates that you should read the user manual and safeguard it.

This symbol identifies information from the manufacturer (see chapter “Manufacturer information”).

This symbol marks the serial number of the device.

This symbol indicates that the middle (tip) of the output plug is positive (+) and the body (ring) of the output plug is negative (-).

IP 22

Protected against ingress of solid foreign objects with a diameter greater than 0.1" (2.5 mm). Protected against access with a finger. Second digit: Protection against dripping water falling vertically when the housing is tilted up to 15°.

Keep out of reach of children.

Hotline (After sales support)

E-Mail (After sales support)

Safety

Explanation of notes

The following symbols and the signal words are used in this user manual.

This signal symbol/word designates a hazard with moderate degree of risk which may lead to death or severe injury if not avoided.

CAUTION!

This signal symbol/word designates a hazard with low risk that, if not avoided, may result in minor or moderate injury.

NOTICE!

This signal word warns against potential damages to property.

General safety information

WARNING!

Danger for children and persons with impaired physical, sensory or mental capacities (e.g. partially disabled persons, older persons with reduced physical and mental capacities) or lack of experience and knowledge (e.g. older children).

- This blood pressure monitor may be used by persons with impaired physical, sensory or mental capacities or those lacking experience and knowledge if they are supervised or have been instructed in how to safely use the blood pressure monitor and have understood

the risks associated with operating it. Children must not play with the blood pressure monitor. Cleaning and user maintenance must not be performed by unsupervised children.

- Keep children and young persons under 18 years of age away from the blood pressure monitor.
- Do not let children play with the packaging. Children could suffocate themselves.

WARNING!

Risk of injury!

Improper use of the blood pressure monitor may result in injury.

- Only use the blood pressure monitor for its intended purpose. Do not misuse in any way.
- To avoid inadvertent strangulation, keep this blood pressure monitor away from children and do not put hoses around the throat.

- This blood pressure monitor cannot serve as a substitute for medical advice or treatment! The measurements taken are only for reference. If you experience problems with your health, it is imperative that you consult your physician! Do not take any medications based on measurements or do not stop taking any medication without having first consulted your physician!
- The blood pressure monitor is not suitable for infants, children or individuals who cannot express their intentions.
- The blood pressure monitor is not suitable for people with severe cardiac arrhythmia.
- Persons with serious circulation problems may experience discomfort. Consult your physician before use.
- Contact your physician if test results regularly indicate abnormal readings. Do not attempt to self-treat these symptoms without consulting a physician.
- Before use, consult a physician if you have any doubts about using the blood pressure monitor.

-
- If you experience any anomalies while using the blood pressure monitor (e.g. allergies), immediately stop using it and consult your physician.
 - Measuring your blood pressure and pulse rate with the blood pressure monitor is not a substitute for medical diagnosis and treatment.
 - Keep all liquids away during use and in particular never immerse the blood pressure monitor in liquids. Otherwise, injury and detrimental health effects may result due to a short-circuit.
 - Do not use the blood pressure monitor in bed or when sleeping.
 - Please consult your doctor before using the blood pressure monitor for the first time if you have any of the following:
 - a.** injuries on your arm, **b.** wounds on your arm that have not yet healed, **c.** illnesses in the arm area, **d.** irregularities in the arm area, **e.** pain with unknown cause in the arm area, **f.** lack of pain perception in the arm area, **g.** metals or implants in the arm area.
-

- For persons with irregular or unstable circulation due to diabetes, liver diseases, arteriosclerosis or other diseases, the blood pressure levels on the wrist may differ from the ones on the upper arm. Nevertheless, monitoring the blood pressure development on the arm or wrist is helpful and important.
- People with vascular constrictions, liver diseases or diabetes, people with pacemakers or weak pulse as well as pregnant women should consult their physician before measuring their blood pressure themselves. Different levels may be produced due to their condition.
- Persons suffering from arrhythmia like atrial or ventricular premature heart beats or atrial fibrillation may only use this blood pressure monitor after consultation with their physician. The oscillometric measurement method can lead to incorrect measurements under certain circumstances.

-
- Excessive measurements may result in injuries to the patient due to blood flow disruptions.
 - The cuff should not be placed on a wound as this may cause other injuries.
 - Do not attach the cuff to a limb that is used for IV infusions or other intravascular accesses, therapies or arteriovenous (A-V) shunts. Inflating the cuff may temporarily block the blood flow and potentially injure the patient.
 - The cuff should not be placed on the arm on the side of a mastectomy (removal of the mammary gland). In case of double mastectomy, use the side of the less dominate arm.
 - The pressure load of the cuff may temporarily cause loss of function for monitoring instruments which are simultaneously used on the same arm.
 - A compressed or kinked connecting hose may cause a continual cuff pressure leading to circulation problems and possibly harmful injuries.

- Ensure that the use of the blood pressure monitor doesn't lead to a permanent impairment of the circulation.
- Prolonged hyperinflation of the cuff may cause an ecchymoma (patchy bleeding) on the arm.
- Only use the approved cuff for this blood pressure monitor. Using other cuffs may lead to incorrect measurements.
- If the system is not stored or used within the temperature and humidity ranges specified by the manufacturer, it may provide incorrect values.
- Avoid excessive pressure when using the cuff to prevent insufficient circulation and blood arrest. Immediately switch the blood pressure monitor off if air is not released from the cuff directly after taking a measurement. In the event of prolonged high pressure in the cuff, the bloods vessels and nerves are constricted, which may result in irreparable damage to tissue and loss of limbs. Avoid falling asleep while taking a measurement.

- Check the blood pressure monitor for damage before each use. If an error or defect (e.g. caused by children, pets or pests) is at hand, the blood pressure monitor must not be used. Please contact the service center or the distributor.
- When performing frequent measurements, some users may suffer from irritation to the skin caused by the cuff during the course of the day. However, this irritation typically disappears by itself after removing the cuff.
- Do not use the blood pressure monitor if it is damaged. Using a damaged device can result in injury, pose serious risks and result in inaccurate measurements.

NOTICE!

Risk of damage!

Improper use of the blood pressure monitor may result in inaccurate measurements and damage to the blood pressure monitor.

- Do not disassemble or attempt to repair the blood pressure monitor.

- Do not insert foreign objects into openings.
- Do not crush the pressure cuff.
- Do not subject the blood pressure monitor to strong magnetic fields (e.g. transformers). Otherwise this could cause errors in readings or damage to the blood pressure monitor.
- Do not use the blood pressure monitor within 4.92 feet (1.5 meters) of a shortwave or microwave appliance or a high-frequency electrosurgical device as this poses a risk of irritation to the skin.
- Do not use the blood pressure monitor in the mountains at altitudes of over 9,843 feet (3,000 meters).
- Please note that portable and mobile HF (high-frequency) communication equipment (e.g. mobile phones) could interfere with medical electrical devices.
- Medical electrical devices are subject to special precautionary measures with respect to the EMC (electromagnetic compatibility). Please therefore comply

with the EMC information provided (see chapter “Information about electromagnetic immunity”) for installing and commissioning the blood pressure monitor.

- Only use the blood pressure monitor at an ambient temperature of 50 F to 104 F (+10 °C to +40 °C) and a humidity of 15 % to 93 %.
- If the blood pressure monitor was stored at temperatures below 32 °F (0 °C), let it rest in a warm location for approx. 15 minutes before use.
- If the blood pressure monitor was stored at temperatures over 104 °F (40 °C), let it rest in a cool location for approx. 15 minutes before use.
- Do not open the blood pressure monitor.
- Protect the blood pressure monitor against dirt, direct sunlight, severe temperature fluctuations, chemicals and mechanical shocks.
- Do not let the blood pressure monitor fall, do not shake it and do not expose it to shocks.

- Only use enclosed or original replacement cuffs to ensure that the measurements are correct.
- Do not start the measurement process until the cuff has been fitted.
- Do not immerse the blood pressure monitor in water and do not rinse it off under running water.
- Protect the cuff and its hose from sharp and pointed objects.
- Do not kink the cuff hose.
- Do not place any heavy objects or objects with sharp edges on the blood pressure monitor, the cuff or cuff hose.
- Do not use the blood pressure monitor at the same time as other medical and electrical devices of any kind.
- Do not use the blood pressure monitor while showering, swimming, visiting the sauna, bathing or when located in an environment with a high degree of air humidity.

Battery notes

WARNING!

Health hazard!

There is a high health hazard if children or animals swallow the batteries.

- Store the batteries so they are not accessible to children and animals. If you suspect a battery has been swallowed, seek immediate medical advice.

WARNING!

Risk of explosion and fire!

Using the blood pressure monitor or the batteries improperly may result in fire or an explosion.

- Do not throw batteries in an open fire.
- Do not charge any non-rechargeable batteries.
- Never destroy the batteries.
- Never short circuit the batteries.

- Promptly remove empty batteries from the battery compartment.
- Only use batteries of the same type (see chapter “Technical data”).
- Do not use old and new batteries together.
- Before inserting the batteries, make sure the contacts inside the battery compartment and on the battery itself are clean and intact. Clean them if necessary.
- Remove the batteries from the blood pressure monitor if it is not being used for more than 3 months to prevent damage caused by leaky batteries.

**WARNING!****Risk of poisoning and chemical burns!**

Ingesting batteries or coming into contact with battery acid can result in poisoning and chemical burns.

- Store the batteries so that they are not accessible for children.
- Batteries may be installed or replaced only by an adult.
- If it is no longer possible to securely close the battery compartment, dispose of the blood pressure monitor as described in the chapter “Disposal”.
- If you suspect that a battery has been swallowed or otherwise incorporated, promptly consult a physician.
- Do not touch leaky batteries. If you do happen to come into contact with battery acid, wash the affected area thoroughly with plenty of clean water.

Safe use of the power supply unit

Risk of electric shock!


Improper use of the power supply unit may result in injury.

- Only use a medical power supply unit which is suitable for the blood pressure monitor (DC 6.0 V, 600 mA, not provided).
- Only use the power supply unit in dry and indoor areas. Do not let the power supply unit come into contact with water or other liquids.
- Connect the power supply unit only to an alternating current (AC) power supply. Make sure that the voltage corresponds with the voltage details on the power supply unit.
- Only connect the power supply unit to a properly installed and easily accessible socket so that you can quickly pull the plug in an emergency.

- Never cover the power supply unit. Ensure that sufficient ventilation is provided in order to prevent the power supply unit from overheating.
- Do not use the power supply unit if the housing or cord is defective or has other visible damages.
- In the event of damage, only have the power supply unit serviced by the manufacturer, its customer service team or a similarly qualified person to avoid risks.
- Always disconnect the power supply unit from the socket if you would like to clean the blood pressure monitor or if you do not use the blood pressure monitor for a prolonged period.


Product description

The blood pressure monitor is only designed for non-invasive measurement of systolic, diastolic blood pressure and pulse rate of adults, using the oscillometric method. The results are shown on the display **5** of the blood pressure monitor **1**. The blood pressure monitor can run on batteries **2** or a power supply unit (not included in the package contents) connected to the power plug jack **12**.

- 4 Cuff
- 5 Display
- 6 M button
- 7 START/STOP button
- 8 SET button
- 9 Battery compartment cover
- 10 Battery compartment
- 11 Air plug jack
- 12 Power plug jack
- 13 Cuff hose
- 14 Air plug

Description of display

- 15 Indicator for systolic pressure
- 16 Indicator for diastolic pressure
- 17 Pulse rate indicator
- 18 Pulse rate/memory slot
- 19 Memory slot indicator
- 20 Date
- 21 Time

- 22 Arrhythmia indicator
- 23 Average indicator
- 24 Air release indicator
- 25 Indicator for measurement classification according to the WHO
- 26 Low battery indicator
- 27 User icon (♂ or ♀)

First use

Checking the blood pressure monitor and package contents

NOTICE!

Risk of damage!

If you are not cautious when opening the packaging with a sharp knife or other pointed objects, you may quickly damage the blood pressure monitor.

- Be very careful when opening it.
- Before first use, use your fingernail to carefully remove the protective foil from the display.

Operation

1. Take the blood pressure monitor out of the packaging.
2. Check to make sure the package contents are complete (see **Fig. on page 2**).
3. Check whether the blood pressure monitor or individual parts are damaged. If this is the case, do not use the blood pressure monitor. Contact the warranty provider as specified on the warranty card.

Operation


Inserting or changing batteries

NOTICE!

Risk of damage!

Incorrect battery installation may damage the blood pressure monitor.

- Make sure that the polarity of the batteries (+ or -) is correct.

1. Open the battery compartment **10** by removing the battery compartment cover **9** (see **Fig. 1**).
2. To replace the batteries **2**, remove the old batteries first. Insert 4 batteries (size AA) in the battery compartment.
3. Close the battery compartment by replacing the battery compartment cover.

Using a power supply unit

The blood pressure monitor can also be operated with a power supply unit (not included in the package contents). The use of a medical power supply unit (DC 6.0 V, 600 mA) is prescribed. The adapter connection pin must be positive inside and negative outside with a 0.08" (2.1 mm) coaxial joint.

1. Connect the output plug of the power supply unit (not included in the package contents) to the power plug jack **12**.
2. Connect the power supply unit to a socket.

System settings

The system settings allow you to select users, and set date, time and time format.

1. Push the SET button **8**.
A user icon ( or  **27**) will show on the display **5**.
2. Push the M button **6** to select either user icon 1 or user icon 2. This allows up to two users to save and accumulate the test results (60 memory slots per user).
3. Push the SET button to confirm the setting. A message “**Yr**” (year) will show on the display.
4. Push the M button to set the year and then push the SET button to confirm.
5. Repeat steps 3 to 4 to set the month, date, hour time and minute time.
The date **20** and time **21** are shown at the lower left corner of the display.
6. After the minute time is set, push the M button to select 12-hour time format or 24-hour time format.

The 12-hour time format is indicated by the message “12” and the 24-hour time format by “24”.

7. Push the SET button to access voice output mode. With the M button, you can configure voice output “L1” (English), “L2” (Spanish) or “OFF” for no voice output.
8. Push the SET button to access volume mode. Set the volume of the voice output with the M button.
9. Push the START/STOP button **7** to save the settings. The blood pressure monitor will also automatically save all settings and switch itself off without any operation for 3 minutes.

It is important to set the date and time. Otherwise, you will not be able to save the measurement results correctly and to access them again later.

Notes on use

For ideal measurement results, follow the basic rules below when measuring your blood pressure.

- Use of the cuff does not pose any biological hazard/risk to the user if it is worn on the skin.
- Before using the blood pressure monitor, let it reach room or ambient temperature to prevent falsified measurement results.
- If possible, always measure your blood pressure at the same time of day for comparable results.

- Test comparisons should only be performed if the measurement is taken on the same arm, in the same position and at the same time of the day.
- Morning time or any time you feel relaxed and stress-free is the best for measurements.
- Always measure blood pressure in a calm and relaxed state, not right after a large meal.
- Do not measure your blood pressure if you are hot or cold, only do so when your body temperature is normal if possible.
- Do not stand during measurements. Sit or lie down in a relaxed position and keep your arm level with your heart (see **Fig. 3 on page 36**).
- Do not move and do not speak while performing a measurement.
- Avoid eating, working out or taking a bath 30 minutes before the test.
- Give the circulatory system about five minutes to calm down each time before a measurement is made.
- If the measurements are not plausible, repeat the measurement.
- Wait at least five minutes before taking another measurement for the same person.
- The results of the blood pressure measurements are for information only and are no substitute for a medical examination. Do not make any medical decisions based on your own measurements (e.g. with respect to taking medication or their dosage). Consult your physician if you have questions about the measurements.
- Do not use this blood pressure monitor to perform any measurements on preeclampsia patients or newborn babies.

-
- Consult a physician before you use the blood pressure monitor on pregnant persons.
 - Only use the blood pressure monitor for persons with an upper arm circumference between 8.66" and 16.54" (22 and 42 cm) (adults).
 - The function of the affected body part may be impaired when inflating the cuff.
 - The cuff hose must not be kinked, compressed or otherwise restricted.
 - Always measure your blood pressure on the same arm to the extent possible – ideally, the left arm. Measuring on the left arm may provide more accurate results as it is closer to the heart.
 - If air escapes from the cuff, replace the cuff with a new one. It is generally recommended to have the cuff replaced in due time in order to ensure correct function and accuracy. Please consult your local dealer or the customer service.
 - Do not use the cuff on limbs whose arteries or veins are undergoing medical treatment, e.g. in the form of intravascular access, intravascular therapy or arteriovenous (A-V) shunt.
 - Do not fit the cuff on persons, who have undergone a mastectomy.
 - Please keep in mind that the restriction of circulation through the inflation of the cuff may temporarily cause other measurement devices connected to the same body part to temporarily malfunction.
 - Avoid performing frequent measurements and prolonged pressure in the cuff to prevent injury.
 - The following illnesses or symptoms may result in incorrect measurements or affect the measurement accuracy of the blood pressure monitor:
-

Performing a measurement

- Diseases of the cardiovascular system,
 - very low blood pressure,
 - diabetes,
 - circulatory disorder and cardiac arrhythmia,
 - shivering,
 - shaking.
- The following conditions may result in incorrect measurements or affect the measurement accuracy of the blood pressure monitor:
- improper positioning of the cuff,
 - improper body position,
 - Feeling anxious or nervous.

Performing a measurement

CAUTION!

Risk of injury!


Improper use of the blood pressure monitor may result in injury.

- Avoid prolonged exposure to the pressure exerted by the cuff, e.g. by kinking the cuff hose and by performing prolonged or frequent measurements.
- Do not place the cuff over wounds to prevent further injury.

- Remove the cuff from your arm in the event of a malfunction of the blood pressure monitor.

Fitting the cuff

1. Expose the arm where you would like to measure the blood pressure. Use your left arm if possible.
2. Wrap the cuff **4** around the exposed upper arm. The bottom edge of the cuff must be $0.39\text{''}-0.79\text{''}$ ($1\text{-}2\text{ cm}$) above the elbow and over the artery (see **Fig. 2**). The cuff hose **13** must point towards the middle of your palm.


3. Now place the free end of the cuff around your arm. The cuff must be fitted tightly enough so that two fingers fit under the cuff.
4. Once the cuff has been properly fitted, close the hook and loop fastener on the cuff.
5. Connect the air plug **14** of the cuff to the air plug jack **11**.


Assuming the correct posture

You can perform the measurement while sitting or lying down. Be absolutely sure to observe the following instructions:

- The cuff must be roughly level with your heart, but no higher (see **Fig. 3**).

- Rest your arms in a relaxed position. A chair with armrests or a table is required for measurements performed while sitting (see **Fig. 4**).

- Do not move and do not speak while performing a measurement.

If you would like to measure your blood pressure while sitting, you must also observe these instructions:

- Sit in a relaxed position when measuring your blood pressure and lean back.
- Place your feet flat on the ground.
- Do not cross your legs.

Starting a blood pressure measurement

- You can always cancel blood pressure measurements that are in progress (e.g. when the pressure in the cuff becomes too extreme) by pushing the START/STOP button.
- You may skip steps 1–3 if you do not intend to save the measurement results into memory.

1. Push the SET button **8**.
A user icon (**1** or **2**) **27** will show on the display **5**.
2. Push the M button **6** to select either user icon 1 or user icon 2.
This allows the measurements to be saved into the memory slots of the selected user after the measurement is done.
3. Push the START/STOP button **7** to save the setting.
4. Push the START/STOP button to start the measurement.
5. You will hear a beep and the advice that the measurement will now begin.
If there is residual air in the cuff **4** from the previous measurement, air will be released from the cuff until the pressure is stabilized before the measurement. During this time, the air release indicator **24** will flash on the display.
6. The cuff will be initially inflated to a pressure of 190 mmHg. You will feel it tightening. If your systolic blood pressure is over 190 mmHg, the blood pressure monitor will automatically re-inflate the cuff to measure the blood pressure.
7. As soon as a pulse is detected, the pulse rate indicator **17** will start to flash. Keep relaxed during this time, avoid speaking and moving any body parts.

8. Once the measurement is complete, the air will be released from the cuff.

The indicator for systolic pressure **15**, diastolic pressure **16** and the pulse rate **18** will be shown on the display. The measurements will be read out.

The corresponding indicator for measurement classification according to the WHO **25** will also be shown (see chapter “WHO classification”).

Deleting or storing measurement results

If the number of measurements exceeds 60 memory slots per group, the most recent measurements will appear at first and the oldest measurements will be deleted as a result.

- If you think the current readings are inaccurate (due to improper conditions or any other reasons), you can delete them by pressing the SET button **8**.
- If the results are not deleted, they will be automatically stored into the memory slots of the selected user (see steps 1–3 in the chapter “Starting a blood pressure measurement”).

Deactivation


- Push the START/STOP button **7** to switch off the blood pressure monitor.
- The blood pressure monitor automatically switches off after 3 minutes.

Evaluating the measurements

About blood pressure


Blood pressure is the force of blood pushing against the walls of arteries. It is typically measured in millimeters of mercury (mmHg).

Systolic blood pressure is the maximum force exerted against blood vessel walls each time the heart beats. Diastolic blood pressure is the force exerted on blood vessels when the heart is resting between beats (see **Fig. 5**).

Example: Fluctuations within a day (male, 35 years of age)

The ABO 523 features an indicator for the blood pressure classification according to the guidelines of the World Health Organization. The following diagram (color coded with respect to the unit of measure) shows the measurement results.

Slightly elevated blood pressure

Evaluating the measurements

An individual's blood pressure frequently changes throughout the course of a day. Excitement and tension can cause blood pressure to rise, while drinking alcohol and bathing can lower blood pressure. Certain hormones like adrenaline (which your body releases under stress) can cause blood vessels to constrict, leading to a rise in blood pressure.

If the blood pressure measuring numbers become too high, it means the heart is working harder than it should.

High blood pressure is a dangerous disease and can affect one's quality of life. It can lead to a multitude of problems such as cardiac insufficiency, kidney failure and cerebral bleeding. A healthy lifestyle and regular visits to one's physician make it much easier to keep high blood pressure and associated diseases under control, particularly if they are diagnosed in the early stages.

If blood pressure is classified as optimal, normal or elevated normal, the green backlight will illuminate. If blood pressure is classified as slightly elevated, the yellow backlight will illuminate.

If blood pressure is classified as moderate high or severe high blood pressure, the red backlight will illuminate.

Cardiac arrhythmia

- If the blood pressure monitor detects an irregular heart rhythm two or more times during the measuring process, the arrhythmia indicator **22** will be shown on the display **5**.
- Irregular heart rhythm is defined as a heart rhythm that is either 25 % slower or faster than the average heart rhythm.
- Consult your physician if the arrhythmia indicator frequently appears with your measurement results. At the same time, make sure you comply with the basic rules

for correct measurement of blood pressure (see chapters “Notes on use” and “Performing a measurement”).

Average of measurement results

- When the blood pressure monitor is switched off, push the M button **6** to show the average values from the last three measurement results of the previously used group.
- The average indicator **23** will show on the display **5** along with the average values.

To check the average values from another user, select the desired user before pressing the M button. See steps 1–3 in the chapter “Starting a blood pressure measurement”.

Checking memory

1. When the blood pressure monitor is switched off, push and hold the M button **6** to access the stored readings for the current user.
The memory slot indicator **19** will be shown along with the current memory slot **18**.
2. Push the M button repeatedly to scroll through all memory slots with the saved measurement results.
There are 60 memory slots for each user. When all memory slots are used, the oldest stored measurement results will be erased to free up a memory slot for the new results.

To check the memory from another user, select the desired user before you push and hold the M button. See steps 1–3 in the chapter “Starting a blood pressure measurement”.

Deleting memory

- When you are checking the stored readings for a user, push and hold the SET button **8** for approximately 3 seconds to delete all stored measurement results.
- You will hear a beep tone when deletion is complete. Memory cannot be recovered once deleted.

WHO classifications

The World Health Organization (WHO) has defined six categories for the classification of blood pressure levels. This classification is provided in the following table:

Range of blood pressure levels	Systoles (mmHg)	Diastoles (mmHg)	Measure
Optimal	< 120	< 80	Check your blood pressure on your own regularly.
Normal	120–129	80–84	
Elevated normal	130–139	85–89	Have your blood pressure checked regularly by a physician.
Slight hypertension	140–159	90–99	
Moderate hypertension	160–179	100–109	Have your physician examine you promptly.
Severe hypertension	≥ 180	≥ 110	

The indicator for measurement classification according to the WHO **25** shows you the range in which the measured blood pressure lies.

If the values for systoles and diastoles are within two different WHO ranges, the device will display the higher of the two ranges.

However, the WHO only offers a general guideline. Variances in the individual blood pressure of different persons are not unusual. Consult your physician to find out what blood pressure levels are normal for you and what ranges are considered dangerous.

Questions and answers

NOTICE!

Risk of damage!

Failure to handle the blood pressure monitor properly could damage it.

- Do not open the blood pressure monitor. Liability and warranty claims are waived in the event of repairs performed by the user or incorrect operation.

Some problems may be caused by minor faults, which you can rectify yourself. To do so, follow the instructions in the following table. If this does not correct the problem, contact the after sales support. Do not repair the blood pressure monitor yourself.

Questions and answers

Question	Answer
What is the difference between a blood pressure measurement taken at home and one at a healthcare facility?	Today, it is assumed that blood pressure measurements at home are more accurate because they better reflect one's everyday situation. Measurements taken at a healthcare facility may be higher. This is referred to as white coat hypertension and can be caused by anxiety or nervousness.
What is the best time of day to measure blood pressure?	In the morning or whenever you feel relaxed and stress-free.
What could be the cause of abnormal measurements?	<ul style="list-style-type: none">- Improper fitting of the cuff.- Make sure that the cuff is fitted neither too tight nor too loose.- Make sure that the cuff is fitted approx. 0.39"-0.79" (1-2 cm) above the elbow.- Improper posture.- Make sure that you keep your body in an upright position.- Anxiety or nervousness.- Take 2 to 3 deep breaths, wait a few minutes and then continue taking the measurement.
What causes different measurements?	Blood pressure fluctuates during the course of the day. Numerous factors can have an effect on blood pressure such as diet, stress, how the cuff is fitted etc.

Question	Answer
Should I put the cuff on my left or right arm? What is the difference?	Blood pressure can be measured on both arms. However, to obtain comparable results, you should always use the same arm. Measurements on the left arm can provide more accurate results as they are closer to the heart.

Troubleshooting

Problem	Possible cause	Solution
Blood pressure results are not within typical range.	The cuff is too tight or not properly positioned on the arm.	Reposition the cuff (approximately 0.39–0.78" (1–2 cm)) above the elbow joint.
	Inaccurate measurement results due to body or blood pressure monitor movement.	Sit in a relaxed position with the arm placed near the heart. Do not move and do not speak while performing a measurement. Make sure that the blood pressure monitor is placed in a stationary position throughout the entire measurement.

FCC Compliance Statement

Problem	Possible cause	Solution
"Err" is displayed.	The cuff fails to inflate properly.	Make sure that the cuff hose is properly connected to the cuff and blood pressure monitor.
	Improper operation	Read the user manual carefully and perform the measurement again properly.
	Pressurization is over 300 mmHg.	
The blood pressure monitor functions unexpectedly during operation.	The blood pressure monitor is defective.	Report this irregularity to the service centre or distributor.

FCC Compliance Statement

Responsible party

Company name: Greater Concepts by Design

Address: 27825 Elie Oak Lane
Spring Texas 77386, USA

Phone number: 1-815-342-5850

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to the product not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Information about electromagnetic immunity

The ABO 523 model satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The ABO 523 model is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the user manual. Portable and mobile HF communications equipment can affect the device.

Information about electromagnetic immunity

Use of the unit in conjunction with non-approved accessories can affect the ABO 523 model negatively and alter the electromagnetic compatibility. The ABO 523 model should not be used directly adjacent to or between other electrical equipment.

Instructions and declaration of manufacturer – electromagnetic emissions

The ABO 523 model is intended for use in the electromagnetic environment specified below. The customer or the user of the ABO 523 model should ensure that it is used in such an environment.

Emission test	Conformity	Electromagnetic environment – instructions
Radiated emission CISPR 11	Group 1, class B	The ABO 523 model only uses RF energy only for its internal function. Therefore, its emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted emission CISPR 11	Group 1, class B	The ABO 523 model is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Instructions and declaration of manufacturer – electromagnetic immunity			
<p>The ABO 523 model is intended for use in the electromagnetic environment specified below. The customer or the user of the ABO 523 model should ensure that it is used in such an environment.</p>			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic discharge IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/burst IEC 61000-4-4	± 2 kV, 100 kHz, For AC power sources	± 2 kV, 100 kHz, For AC power sources	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV (differential mode)	± 0.5 kV, ± 1 kV (differential mode)	Mains power quality should be that of a typical commercial or hospital environment.

Information about electromagnetic immunity

<p>Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11</p>	<p>0 % UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles</p>	<p>0 % UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles</p>	<p>Mains power quality should be that of a typical commercial or hospital environment.</p>
<p>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</p>	<p>30 A/m; 50 Hz or 60 Hz</p>	<p>30 A/m; 50 Hz or 60 Hz</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>

Information about electromagnetic immunity

<p>Conducted RF IEC 61000-4-6</p>	<p>3 V for 0.15–80 MHz; 6 V in ISM and am- ateur fre- quencies between 0.15 and 80 MHz</p>	<p>3 V for 0.15–80 MHz; 6 V in ISM and am- ateur fre- quencies between 0.15 and 80 MHz</p>	<p>Portable and mobile RF communications equip- ment should be used no closer to any part of the ABO 523 model, in- cluding cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p>
<p>Radiated RF IEC 61000-4-3</p>	<p>385 kHz, 27 V/m 450 MHz, 28 V/m 710 MHz, 745 MHz, 780 MHz 9 V/m 810 MHz, 870 MHz, 930 MHz 28 V/m 1720 MHz, 1845 MHz, 1970 MHz 28 V/m 2450 MHz, 28 V/m 5240 MHz, 5500 MHz, 5785 MHz 9 V/m</p>	<p>385 kHz, 27 V/m 450 MHz, 28 V/m 710 MHz, 745 MHz, 780 MHz 9 V/m 810 MHz, 870 MHz, 930 MHz 28V/m 1720 MHz, 1845 MHz, 1970 MHz 28 V/m 2450 MHz, 28 V/m 5240 MHz, 5500 MHz, 5785 MHz 9 V/m</p>	<p>Recommended separation distance</p> $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ <p>800 MHz to 2.7 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range.</p> <p>The product may experi- ence interference when it is near devices with the following symbol:</p> 

Information about electromagnetic immunity

Recommended separation distances between portable and mobile RF communications equipment and the device

The ABO 523 model is intended for use in an electromagnetic environment in which radiated disturbances are controlled. The customer or the user of the ABO 523 model can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ABO 523 model as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1}\right]\sqrt{P}$	800 MHz to 2.7 GHz $d = \left[\frac{7}{E_1}\right]\sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTICE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTICE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Cleaning and maintenance

NOTICE!

Risk of damage!

Improper handling of the blood pressure monitor may result in damage to the blood pressure monitor.

- Do not use any aggressive cleaners, brushes with metal or nylon bristles or sharp or metallic cleaning utensils such as metal sponges and the like. They could damage the surfaces.
- Do not in any case place the blood pressure monitor in the dishwasher. It would be destroyed as a result.
- Switch the blood pressure monitor off when cleaning it.
- Do not submerge the blood pressure monitor in water or other liquids in order to clean it.
- Do not soak the cuff in water.

Cleaning

The blood pressure monitor and cuff must be cleaned before being passed on to other persons to prevent the transmission of diseases.

1. Always remove the batteries **2** from the blood pressure monitor before cleaning it.
2. Carefully clean the surfaces of the blood pressure monitor and the cuff **4** at regular intervals (e.g. after every 20th use) with a soft, damp cloth and a small amount of rubbing alcohol. Please ensure that no moisture enters the device. A mild cleaning agent may be used for stubborn dirt.
3. After cleaning, use a soft cloth to dry all parts completely.
4. Let the cuff air dry at room temperature.

Cleaning the cuff

- Clean the surface of the cuff **4** with a soft cloth and a small amount of rubbing alcohol.
- Wipe the cuff off with a damp cloth (water).
- Let the cuff dry at room temperature.

Changing the batteries

When the low battery indicator **26** shows, change the batteries **2** according to the steps described in the chapter “Inserting or changing batteries”.

Storage

If you do not use the blood pressure monitor for a prolonged period of time, please observe the following notes:

- The blood pressure monitor doesn't require any maintenance.

- If you do not intend to use the blood pressure monitor for more than 3 months, remove the batteries **2** from the battery compartment **10** or disconnect the device from the power supply unit.
- Store the blood pressure monitor in the enclosed storage bag **3** to protect it.
- Do not store the blood pressure monitor near heat sources.
- Keep the blood pressure monitor away from children, pets and pests.
- Do not expose the blood pressure monitor to direct sunlight, fire, dirt, lint, dust, water, heat and extreme temperatures to prevent damage to the blood pressure monitor.
- Observe the information concerning the storage environment (see chapter "Technical data").
- Do not place any heavy objects on the blood pressure monitor, the cuff **4** or cuff hose **13**.
- Do not allow the blood pressure monitor to fall, do not expose it to shocks and do not throw it.
- Expected service life: Approx. 3 years with 10 measurements per day.

Technical data

Product description:	Blood pressure monitor
Model:	ABO 523 (fully automatic)
Article number:	808487
Display:	Digital LC display, tricolor backlight
Size of display:	4.02" × 2.71" (102.1 × 68.9 mm)
Measurement method:	Oscillometric method
Application part:	The cuff is the application part.
Measurement range:	Systolic blood pressure 60 mmHg–280 mmHg Diastolic blood pressure 30 mmHg–200 mmHg Blood pressure 0 mmHg– 300 mmHg Pulse 30–180 beats/minute
Measurement accuracy:	Pressure: ± 3 mmHg, ± 2 % over 200 mmHg Pulse: ± 5 %
Pressurization:	Automatic
Memory:	120 memory slots in two groups with date and time
Function:	Irregular heartbeat detection WHO classification indicator Average of last 3 measurements Low battery detection Automatic deactivation

Power supply:	4× AA batteries or medical power supply unit (DC 6.0 V, 600 mA, not provided)
Battery life:	approx. 2 months with 3 measurements per day
Service life of the device:	approx. 3 years with 10 measurements per day
Weight (device):	approx. 0.25 lb (357 g) (without batteries)
Weight (cuff):	approx. 0.27 lb (123 g)
Dimensions (L × W × H):	approx. 2.6 × 2.8 × 1.2" (164 × 111 × 60 mm)
Universal cuff (W × L):	approx. 5.71" × 22.05" (145 × 560 mm)

LOT

Batch code, V3821AB0523

SN

Serial number, SN 00001 (sequential number)

2021-08

Date of manufacture, 2021-08 (year, month)

Operating temperature: 50 °F–104 °F (10 °C–40 °C)

Operating humidity: 15 %–93 % R.H.


Working pressure: Atmospheric pressure

Storage temperature: -13 °F–158 °F (-25 °C–70 °C)

Humidity during storage and transport: 15 %–93 % R.H.

Atmospheric pressure during storage and transport: 700–1,060 hPa

Ingress protection rating:

IP 22

Meaning of the first digit: Protected against ingress of solid foreign objects with a diameter greater than 0.1" (2.5 mm). Protected against access with a finger. Second digit: Protection against dripping water falling vertically when the housing is tilted up to 15°.

Classification:

Device with internal power supply, type BF ,

Type BF (body float) protection against electrical shock.

Type BF application equipment with greater protection against electric shock to the body, but not directly to the heart.

This blood pressure monitor complies with the EU regulations and bears the CE mark "CE 0197". This blood pressure monitor complies with the following guidelines, laws and standards (included but not limited):

Safety standard:

EN 60601-1 Medical electrical equipment part 1: General requirements for safety

EMC standard:

EN 60601-1-2 Medical electrical equipment part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility – requirements and tests

Full display reference

Fig. B shows the display **5** with all its indicator segments. You can use this figure to determine whether the display is properly rendering the measurements.

Manufacturer information

**Manufacturer:**

JOYTECH Healthcare Co., Ltd.
No. 365, Wuzhou Road
Yuhang Economic Development Zone
Hangzhou City
311100 Zhejiang
P.R. China

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Disposal

Disposing of the packaging

Dispose of the packaging separated into single type materials. Dispose of paperboard and cardboard with wastepaper and plastics with recyclable waste.

Recycling code PAP 21 means the material is to be classified as “Mixed paper” and is recycled as packaging.

Disposal

Disposing of the packaging

- Dispose of the packaging separated into single type materials. Dispose of paperboard and cardboard with wastepaper and plastics with recyclable waste

Disposing of the blood pressure monitor

Old devices must not be disposed of with household waste!

- Should the blood pressure monitor no longer be capable of being used at some point in time, **dispose of it in accordance with the regulations in force in your federal state or country.** This ensures that old devices are recycled in a professional manner and also rules out negative consequences for the environment.

- Please ensure your recycling information applies to the local regulations and EPA recommendations (www.epa.gov).

Batteries and rechargeable batteries do not belong in the household waste!

As the end user, you are required by law to bring all batteries and rechargeable batteries, regardless of whether they contain harmful substances* or not, to a collection point run by the communal authority or borough or to a retailer so that they can be disposed of in an environmentally friendly manner.

*labeled with: Cd = cadmium, Hg = mercury, Pb = lead

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