User Manual

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
TENS, EMS AND MESSAGE DEVICE

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




Contents

Overview	3
Use	4
Package contents/device parts	6
General information	7
Reading and storing the user manual.....	7
Explanation of symbols	7
Basic information.....	9
Safety	10
First steps	24
Initial inspection of package	24
Adhering the adhesive electrodes to the skin	25
Changing the battery and battery information.....	27
Overview of device functions	29
Application	31
Program overview.....	31
Switching on the device	31
Procedure for fixed programs	
TENS 1–27, EMS 1–7 and Massage 1–10.....	32
Procedure for fixed programs	
TENS 28–30 and EMS 8–10	34
Steps for using the quick start programs.....	35
Preferred Program	36
Warnings, alarms and misapplications	37
Application programs	38
TENS application programs.....	38
EMS application programs	41
Massage program.....	43
Adhesive electrodes	44
Instructions for applying the adhesive electrodes	44
Placement of electrodes.....	45
Placement of electrodes for EMS	
(electrical muscle stimulation).....	47
Technical malfunctions	52
Electrical interference immunity	54
Cleaning and care of the TENS/EMS device	60
Notes on using the TENS/EMS device	61
Technical data, symbols, pictograms	62
Disposal	64

Package contents/device parts

- 1 TEN 603 EMS/TENS device
- 2 Belt clip
- 3 3× AAA batteries
- 4 4× adhesive electrodes
- 5 2× lead wires
- 6  button
- 7  button
- 8 CH2 intensity button
- 9  button
- 10  button
- 11  button
- 12  button
- 13 CH1 intensity button
- 14  button
- 15 Channel CH1 output jack
- 16 Channel CH2 output jack
- 17 Quick start mode
NECK/BACK/KNEE
- 18 Treatment duration in minutes
- 19 Keypad lock symbol
- 20 Low battery symbol:
replace batteries
- 21 Output intensity for
channel CH2 (right)
- 22 Pulse length displayed in μ s
- 23 Pulse frequency displayed
in Hz
- 24 Output intensity for
channel CH1 (left)
- 25 SET: Device is in
programming mode
- 26 Current program number
- 27 Therapy mode:
TENS/MASSAGE/EMS
- 28 Battery compartment cover
- 29 User manual

General information

Reading and storing the user manual

This user manual accompanies this TENS, EMS and massage device (hereinafter referred to only as “TENS/EMS device”). It contains important information on start-up and handling.

Before using the TENS/EMS device, read the user manual carefully. This particularly applies to the safety instructions. Failure to heed this user manual may result in severe injury or damage to the TENS/EMS device.

Store the user manual for future use. Make sure to include this user manual when passing the TENS/EMS device on to third parties.

Explanation of symbols

The following symbols and signal words are used in this user manual, on the TENS/EMS device or on the packaging. The safety symbols displayed in this user manual are instructions on how to properly use the TENS/EMS device and are intended for your safety.

This signal symbol/word indicates a hazard with high risk that, if not avoided, results in death or a severe injury.

WARNING!

This signal symbol/word designates a hazard with moderate risk that, if not avoided, may result in death or severe injury.

CAUTION!

This signal symbol/word designates a hazard with low risk that, if not avoided, may result in minor or moderate injury.

Read and observe the user manual!

Warning/danger: Improper use carries the risk of serious injury, damage and risk to life!

These instructions must be followed!

Warning/danger: The device must not be used by persons with pacemakers!

IP 22 The first number 2: represents protection against hard foreign objects with a diameter of .5 inches (12.5 mm) or larger. The second number 2: represents protection against drops of water falling vertically when the housing is tilted by up to 15 degrees. Drops of water falling at a perpendicular angle have no effect when the housing is tilted by up to 15 degrees.

Basic information

What are the features of a TENS/EMS device?

The TEN 603 TENS/EMS device is a device for electrical stimulation. It conducts an electric current to the skin.

EMS (electrical muscle stimulation) = Muscle tissue is stimulated electrically.

TENS (transcutaneous electrical nerve stimulation) = The nerves are stimulated electrically.

This handy, convenient TENS/EMS device is used to electrically stimulate targeted areas of the body.

Information on using TENS/EMS

The principle purpose of the TENS/EMS device is to stimulate the nerves and muscles.

How does TENS work?

According to scientific theory, electrical stimulation therapy works in several ways:

- The gentle electrical pulses travel through the skin to the nearby nerves, blocking pain or rather preventing it from reaching the brain from its point of origin.
- The gentle electrical pulses increase the body's own production of natural pain-reducing substances such as endorphins.

- Furthermore, it can be assumed that the electrical stimulation also stimulates blood flow. Muscles contract upon contact with electrical stimulation and then expand. The repeated contracting and relaxing causes the blood to flow in and out, improving circulation.

How does EMS work?

The EMS device sends soothing pulses through the skin and stimulates the nerves in the area being treated. When the muscle receives this signal, it contracts as if the brain had sent the signal itself. As the strength of the signal increases, the muscle contracts in the same way as it does during exercise. If the signal is interrupted, the muscle relaxes and the process repeats. The goal of electrical muscle stimulation is contraction or vibration in the muscles. Normal muscle activity is controlled by the central and peripheral nervous system, which sends electrical signals to the muscles. The EMS functions in a similar way, but it uses an external source (stimulator) in combination with electrodes placed on the skin to conduct the electrical pulses into the body. The pulses stimulate the nerves so that they send targeted signals to specific muscles, which then react by contracting – just like during normal muscle activity.

Safety

General safety information

- If there is a defect, you must not repair, modify (change) or use the TENS/EMS device. When used incorrectly, the stimulation current can lead to pain, injury or burns.

- If you experience skin changes, pain or swelling, feel unwell, or experience any other irregularities while using the TENS/EMS device, stop using it immediately and consult your physician. Before using, remove all metal objects such as jewelry, belts, wristwatches or other accessories from your body in order to prevent these from coming into contact with the TENS/EMS device or the adhesive electrodes. Do not use the TENS/EMS device while operating a car; do not undertake any other activities during use of the device.
- If you are not sure how to use the TENS/EMS device, you should consult your physician before using it.
- In cases of unexplained pain, swollen muscles or after a severe muscle injury, only use the TENS/EMS device after consulting your physician. Use of the TENS/EMS device cannot replace medical diagnosis and treatment.
- Keep this user manual for the entire service life of the product for later consultation and include it if you pass the TENS/EMS device on to third parties. Make sure that the user manual is also accessible to third parties. The user manual is a component of the TENS/EMS device.
- Misuse and use other than the intended use must be avoided.

- No equipment from other devices may be used.
- If irregularities arise during use, stop using the device immediately.
- Neither the connecting cords nor the adhesive electrodes may be bent at sharp angles.
- Do not place any heavy or sharp-edged objects on the TENS/EMS device or the adhesive electrodes.
- Check the device and the adhesive electrodes for problems before every use. If a problem or defect is present, the device and adhesive electrodes must not be used.
- Never wear any body jewelry or adhesive tattoos in the area being stimulated during use of the TENS/EMS device.
- Do not use the device on areas of the body that are tattooed.

Proper use

Which use/application is the TENS/EMS device suited for / in which environment should it be used?

- Only use the TENS/EMS device for its intended purpose, which is the external application of low frequencies (electrical pulses in the range of up to 150 Hz and up to 50 V) to the human body.

- The TENS/EMS device is exclusively intended for external use (application on the skin) in people with the intention of stimulating the nerves and muscles.
- The device may be used as part of TENS therapy to treat pain in cases of acute pain, arthritis, rheumatoid pain and other conditions causing chronic pain.
- Unless instructed otherwise by a physician, we recommend an average treatment time of 30 minutes, up to three times a day.
- The sensation of intensity varies widely from day to day and can be adjusted by the user to suit their needs through the TENS/EMS device's intensity control settings.
- TENS pain treatment using the TENS/EMS device is purely symptomatic in nature. It can contribute to a reduction in pain and, under some circumstances, can aid recovery. If you have questions about the therapy, consult your physician.
- EMS applications using the TENS/EMS device work to stimulate the muscles.

! ⚠ Which areas of application and use are the TENS/EMS device not suited for / in which setting should it not be used?

- The TENS/EMS device may not be used simultaneously with other medical and electrical devices of any kind.
- Do not use the TENS/EMS device while showering, swimming, using a sauna, bathing or in any other environment with high humidity. Keep all liquids away during use. Otherwise, injury and detrimental health effects may result due to increased stimulation or a short-circuit. Attention! Risk to life!
- Do not use the TENS/EMS device in bed and/or when sleeping. Do not use the TENS/EMS device if easily flammable substances and gases or explosives are nearby.
- During use, the TENS/EMS device may interfere with other electric or electrical equipment. Therefore, do not use the TENS/EMS device in the vicinity of other electric devices.
- Do not use the TENS/EMS device within 1.5 meters of a shortwave or microwave appliance or a high-frequency electrosurgical device, as this may result irritation or burning to the skin under the electrodes. Do not use the TENS/EMS device at altitudes of over 9,843 feet (3,000 meters).

-
- The TENS/EMS device is meant for personal use only and is not intended for commercial or industrial use.
 - Please note that portable and mobile HF (high-frequency) communication equipment (e.g. mobile phones) may interfere with medical electrical equipment. Medical electrical devices are subject to special precautionary measures with respect to the EMC (electromagnetic compatibility). Therefore, please observe the EMC notes concerning the installation and start-up of the TENS/EMS device.
 -  You must not use the TENS/EMS device under the following circumstances:

- a.** if you have cardiac disease or cardiac arrhythmia (use can lead to cardiac arrest),
- b.** directly on a wound,
- c.** during pregnancy, near the uterus or during labor,
- d.** near the eyes,
- e.** if you have a pacemaker,
- f.** on areas of the body with poor circulation,
- g.** if you have an emotional or mental disorder,
- h.** if you have been diagnosed with dementia (mental deterioration),
- i.** if you have a low IQ (intelligence quotient).

- Under the following circumstances, you must consult your physician before using the TENS/EMS device:
 - a.** acute illness,
 - b.** tumors,
 - c.** infectious disease,
 - d.** fever,
 - e.** blood-pressure conditions,
 - f.** skin conditions,
 - g.** following an accident,
 - h.** nausea or dizziness,
 - i.** disease outbreak,
 - j.** as soon as irregularities arise,
 - k.** pain from an unexplained source,
 - l.** diabetes,
 - m.** seizure disorders,
 - n.** while menstruating,
 - o.** if certain areas of the body are insensitive to pain,
 - p.** if you have metal or implants within your body.
- In cases of non-specific pain such as general headaches, treatment with the TENS/EMS device will be ineffective.
- Do not use the TENS/EMS device if you could injure yourself in any way when startled suddenly.

- The adhesive electrodes of the TENS/EMS device may not be used on open wounds, sensitive areas of the skin, or fresh scars.
- The following persons should not use the TENS/EMS device: children, persons at increased risk, persons with allergies, persons with an immunodeficiency, persons experiencing pain from a unexplained cause, persons with diabetes or circulatory system diseases, persons with poor blood circulation in the exterior arteries and tissues or severe cardiovascular disease. If in doubt, ask your physician.
- The TENS/EMS device must not be used on persons who cannot properly perceive electric stimulation currents. Children are sensitive to stimulation currents! Frail persons and those with disabilities may not be able to notify others when the stimulation current is too intense.

Use by children and young people

- This TENS/EMS device may not be used on children.
- The TENS/EMS device must be kept in a location that is inaccessible to children and young people under the age of 18.
- Ensure that the TENS/EMS device does not fall into the hands of children. Children could ingest the small parts and choke. Children could injure themselves when using the device.

Notes on using the TENS/EMS device

- The adhesive electrodes may only be connected to the TEN 603 TENS/EMS device. Please ensure that the device is always switched off when attaching or removing the adhesive electrodes.
- Always switch off the device before repositioning the adhesive electrodes of the TENS/EMS device during use.
- Under some circumstances, use of the TENS/EMS device can irritate the skin. If skin irritations arise such as redness, blistering or itching, discontinue use of the TENS/EMS device. Do not always attach the adhesive electrodes in the exact same locations, as this can lead to skin irritation. Before using the device, you should thoroughly cleanse and then dry the areas of the skin where you intend to place the adhesive electrodes. The skin should be clean and free of oils. Do not connect the lead wires and adhesive electrodes unless the TENS/EMS device is switched off.
- To avoid damaging the wires, do not pull on them to disconnect them from the adhesive electrode plug connections or from the device. Always pull from the plugs when removing the lead wires from the adhesive electrodes or the device.

-
- The adhesive electrodes may be connected to or disconnected from the device using the sockets and the electric lead wire.
 - Every person responds differently to electrical nerve stimulation. If the treatment is unsuccessful, consult your physician. Please remove the protective film prior to attaching the adhesive electrodes. The adhesive strength of the electrodes depends on the characteristics of the skin, storage, and number of times they have been used. If the full surface of the adhesive electrodes no longer adheres to the surface of the skin, they need to be replaced with new adhesive electrodes. The entire surface of the adhesive electrodes must be in contact with the skin in order to avoid a high local current density, which can lead to burning of the skin. After use, reattach the adhesive electrodes to the protective film and keep them in the polybag to prevent them from drying out. This maintains their ability to adhere.
 - Avoid touching the adhesive electrodes during use, as this can cause a short-circuit under some circumstances, resulting in an excessive current density. This excessive current density can lead to burning and injury!

- DITTMANN textile electrodes, back belt, knee cuff or neck collar, available as accessories, are suitable for use with the TEN 603 TENS/EMS device for the application of special programs for different body areas.
- When using these articles with the TEN 603 TENS/EMS, please observe the user manuals for these articles, particularly their safety instructions.

Where may the adhesive electrodes be applied?

- Every person responds differently to electrical nerve stimulation. The electrodes therefore need not be kept to specific areas. If the application is unsuccessful, consult your physician in order to determine what positioning techniques are best suited for you.
- Familiarize yourself with the proper placement of the adhesive electrodes using the instructions in the chapter “Adhering the adhesive electrodes to the skin” and in the example illustrations in the chapter “Adhesive electrodes” for both TENS and EMS use.
- Do not use adhesive electrodes which are smaller than 1.57×1.57 in or 2.46 in² (40×40 mm, or 16 cm²), as this could result in an excessively high electricity density and cause injury.

-
- The adhesive electrodes may not be modified in size, e.g. through cutting off parts. The distance between electrodes should not be less than 2 in (5 cm) and should not exceed 6 in (15 cm).
 - Do not pull on the wires when removing the adhesive electrodes. Lift the adhesive electrodes by their edges and gently pull off. Ensure when using the TENS application that the electrodes are placed in an area which encloses the pain area. With pain in muscle groups, adhere the electrodes so that they enclose all the relevant muscles.
 - When using an EMS application, please note the following regarding the placement of the electrodes: If you wish to activate surface muscles, adhere the adhesive electrodes parallel to the orientation of the muscle fiber. To reach deep muscle layers, it is recommended that the adhesive electrodes be adhered perpendicular (crosswise) to the orientation of the muscle fiber.

 **In which areas should use of the adhesive electrodes be avoided?**

- The adhesive electrodes may not be used on areas of the body with irritated skin, nor on open and fresh wounds and fresh scars.

Do not adhere the adhesive electrodes to the following areas of the body:

- a. on or in the mouth,
 - b. the eyelids,
 - c. the front of the neck,
 - d. the larynx,
 - e. the throat area,
 - f. the carotid artery,
 - g. near the heart,
 - h. the genitals (sexual organs: penis, testicles, etc.),
 - i. the fingers,
 - j. a pacemaker.
- The adhesive electrodes must not be attached in such a way that the current flows directly through the brain, such as placement on both temples.

The transmission of electric current may lead to lethal heart rhythm disturbances if the current flow is applied to both sides of the thorax (sides, back, or front) or over the chest/heart.

Do not adhere the adhesive electrodes to the head and face (e.g. mouth or eyes).

– Do not use on the neck (especially the sinus caroticus) or the throat, as this may cause severe muscle spasms and closing of the trachea, leading to breathing difficulties or a negative impact on heart rhythm or blood pressure.

Contra-indications

- Avoid treatment in a single area for an extended period of time (more than 30 minutes, up to 3× per day), as the muscles in that area may tire and become sore.
- Skin irritation and burning may occur beneath the stimulation electrodes.
- Electrical stimulation in the vicinity of the eyes or on the head and face may cause headaches and other painful sensations.
- If you experience adverse effects, stop using the device and consult your physician.
- The use of electrodes that are too small or improperly placed may lead to discomfort or burning of the skin.

First steps

The following sections provide precise instructions on the reception, initial inspection, set-up, and operation of the combined stimulation device. For safe and optimal use, we recommend setting up and using the stimulator according to the instructions.

Initial inspection of package

Risk of injury!

Medical electrical devices require special precautionary measures with regard to EMC and must be installed and operated in accordance with the EMC information in the tables at the end of this user manual.

- Remove the device from the packaging and check the contents for damage. If no damage is ascertained, check to verify that the system components are intact and that all components are present in accordance with the user manual. If the contents are damaged or incomplete, inform the shipping carrier and your area retailer.

Adhering the adhesive electrodes to the skin

Risk of injury!

Never apply the electrodes in such a way as to allow electrical current to flow near the heart.

Excessive contact with the adhesive electrodes may cause skin irritation.

- Use only the supplied manufacturer's wires and electrodes.

Risk of injury!

Do not continue to use electrodes if their adhesive strength has weakened. This may cause ineffective treatment or skin irritation.

1. Remove the electrodes **4** from the packaging and peel off the protective film. Store the protective film for later use.
2. Position the electrodes precisely on the correct skin area. (For proper placement, see chapter "Adhesive electrodes".)

Important:

The exact placement on the skin should be determined by a trained physical therapist, or the recommended placements in these instructions should be observed.

The following instructions apply in selecting the distance between electrodes:

- Optimum distance: about 2–6 in (5–15 cm);
 - Less than 2 in (5 cm): Strongly stimulated surface tissue;
 - More than 6 in (15 cm): Large, deep-lying areas receive weak stimulation.
3. Press the edge of the electrode firmly with your finger and lightly press the entire electrode so that full surface contact with the skin is established.
 4. Connect the connecting cables **5** with the electrodes (see Fig. 1).
 5. Hold the wire at the plug and insert into the upper jack **15** / **16** of the device as far as it will go (see Fig. 2).


Fig. 1


Fig. 2

For further information concerning polarity requirements please see chapter “Adhesive electrodes”.

Changing the battery and battery information

DANGER from batteries!

If swallowed, batteries could pose a risk of death.

- Therefore, store the batteries so that they are not accessible to small children.
- If a battery is swallowed, you must seek immediate medical attention.
- If battery acid has leaked, avoid contact with the skin, eyes and mucous membranes. Rinse the affected area immediately with a large amount of clean water and promptly consult a physician or seek medical treatment.
- Batteries must not be recharged (with the exception of rechargeable batteries), taken apart, thrown into fire, or short circuited.
- Protect batteries from excessive heat.
- Remove batteries from the device when they are empty or when the device will not be used for a long period of time. Doing so will prevent damage through leakage.
- Replace batteries regularly.

- Do not mix battery types, brands, rechargeable with non-rechargeable batteries, or batteries with different capacities.
- Insert 3 batteries (AAA), observing the correct polarity (+ and - pole).
Battery types: The TEN 603 TENS/EMS device uses AAA alkaline batteries.
- Do not use any rechargeable batteries!
 1. First remove the belt clip **2** in the direction of the arrow.
 2. To open the cover fastener, press on the arrow-shaped hatched surface of the battery compartment covering **28** and slide outwards to remove it.
 3. Remove the empty batteries **3**.
 4. Insert three new alkaline batteries (AAA).
 5. When inserting the batteries, make sure that their polarity is correct (see marking/stamped information in the battery compartment).
 6. Put the battery compartment covering back on and slide it in the direction of the belt clip holster until it locks back into place.

Overview of device functions


Fig. 3

6 button

- Hold the button for approx. 2 seconds to lock the keypad, in order to prevent unintentional use (with the exception of the  button).

7 button

- Press to select therapy program.
- Hold the  button down for approx. 2 seconds to open the setting mode.

Navigation between treatment time, frequency, pulse width.

13 / 8 CH1-/CH2-intensity button

- Increase ▲ or reduce ▼ the output of the CH1/CH2 channel.
- Increase ▲ or reduce ▼ treatment time, frequency and pulse width for individualized program.

9 button

- Press button to turn on device.
- Press to pause treatment.
- Hold the  button down for approx. 3 seconds to turn off the device.

10 button

- Press this button for a quick start of the knee program.

11 button

- Press this button for a quick start of the lower back program.

12 button

- Press this button for a quick start of the neck program.

14 button

- Press this button to navigate between TENS, MASSAGE and EMS therapy modes.
- Hold the  button down for approx. 2 seconds in order to set the start and finish mode for your personalized program.

Application

Program overview

The TEN 603 TENS/EMS has a total of 68 programs:

- 39 TENS programs (27 fixed programs, 3 adjustable programs, 9 quick start programs – 3 of which are adjustable).
- 19 EMS programs (7 fixed programs, 3 adjustable programs, 9 quick start programs – 3 of which are adjustable).
- 10 massage programs.

The intensity of both channels can be individually selected in all programs.

The following 12 programs can be individually customized:

- **TENS:** Programs P28, P29, P30, Neck Pb3, Back Pb3, Knee Pb3.
- **EMS:** Programs P08, P09 and P10, Neck Pb3, Back Pb3, Knee Pb3.

Switching on the device

- Press the  button  in the middle of the keypad to switch on the device.
After approx. 2 seconds, the device goes into standby mode.
- Please observe the following section for operating the device.

Procedure for fixed programs TENS 1-27, EMS 1-7 and Massage 1-10

1. Connect the electrodes **4** to the device **1** and follow the instructions in the chapter “Adhering the adhesive electrodes to the skin”. Affix the adhesive electrodes to the desired area (see chapter “Adhesive electrodes”).
2. Press the **⏻** button **9** and switch on the device.
3. Press the **Ⓜ** button **14** to navigate between therapy modes EMS, MASSAGE and TENS (Fig. 4, example TENS display)
4. Select the desired program by pressing the **Ⓟ** button **7** (Fig. 5, example of EMS display, Program 02).


Fig. 4


Fig. 5

5. At the start of stimulation, the output intensity is also set by default at 0. Select the desired output intensity with the intensity buttons **8** / **13** for the corresponding channel. For CH1 and CH2, intensity can be individually adjusted (Fig. 6).


Fig. 6

Important:

As soon as you set the intensity setting to 1, that channel is active and current will flow between the adhesive electrodes.

6. We recommend using the locking function, so that buttons are not unintentionally pressed during use. To lock/unlock the keypad, hold the **⏏** button **6** down for approx. 2 seconds (Fig. 7).

Fig. 7

7. After the stimulation time has ended, hold the **⏏** button down for approx. 3 seconds to switch off the device.

Notice:

If you wish to end the stimulation treatment prematurely, press the **⏏** button. The device will shut off automatically if idle for more than 1 minute after being switched on.

Procedure for fixed programs TENS 28–30 and EMS 8–10

TENS 28–30 and EMS 8–10 are preset programs for which you can individually adjust the frequency, pulse width and treatment time.

1. Place the electrodes **4** so that they enclose the pain area, and switch on the device. Select an individual TENS/EMS program (Fig. 8, example TENS program 28) as described in steps 1–4 of the previous chapter.
2. Hold the **P** button **7** down for approx. 2 seconds to activate the setting mode.
3. Set the frequency:

The frequency can only be adjusted for programs TENS 28 and TENS 30. Press the intensity buttons **8** / **13** to set a desired frequency between 1 Hz and 150 Hz (Fig. 9, example program TENS 28).

Confirm your selection with the **P** button and you can then adjust the pulse width (except for in program 29).

Fig. 8

Fig. 9

4. Adjust the pulse width:

You can also select a pulse width from 50 μ s to 250 μ s using the intensity button (Fig. 10, example program TENS 28). Confirm your selection with the P button and continue to select the treatment length.

Fig. 10

5. Set the treatment length:

Press the intensity button to select a treatment duration between 5 and 100 minutes (Fig. 11, example program TENS 28).

Press the P button **9** to confirm your selection and deactivate programming mode.

Fig. 11

6. Start and finish the treatment by following steps 5 through 7 in the previous chapter.

Steps for using the quick start programs

There are 3 quick start programs (Pb1, Pb2 and Pb3) in both TENS and EMS mode for the neck, lower back and knee areas.

The quick start programs are selected with the following buttons:

- **12**  button: Neck
- **11**  button: Lower back
- **10**  button: Knee

1. Position the electrodes **4** on the desired area, for example near the pain area on the neck, and switch on the device. Switch on the device by pressing the P button **9**.

2. Press the **M** button **14** to set TENS or EMS mode (see Fig. 12, example TENS program).
3. Press the **Person** button **12** to select the program Pb1 (see Fig. 12, example TENS program Pb1).

Important: Repeatedly pressing the **Person** button will take you to programs Pb2 and Pb3.

Fig. 12

4. Begin treatment, following steps 5 through 7 in the chapter “Procedure for fixed programs TENS 1–27, EMS 1–7 and Massage 1–10”.

Important:

Frequency, pulse width and treatment time can only be set in the relevant Pb3 programs when steps 3 to 6 in the previous chapter are followed.

Preferred Program

The Preferred Program is a special setting with which you can easily call up a program that has been personalized specifically for you. Your program settings are loaded and activated immediately after the device is switched on. Settings for this individual program can be made according to recommendations from your physical therapist or physician, for example.

Setting your preferred program

1. Select your program and the corresponding settings as described in the previous chapters.

At the start of stimulation treatment, the output intensity for both channels CH1 and CH2 is always set by default at 0.

- Before you set the output intensity, hold the **M** button **14** down for approx. 2 seconds to turn on the preferred program.
The storing of the current program settings in preferred program mode will be confirmed by a long beeping sound.

Deleting the preferred program

- To delete the preferred program, hold the **14** button **M** down for approx. 2 seconds.
- Output intensity must be set to 0 mA before the program can be deleted.
The deletion of the preferred program will be confirmed by a double beeping sound.

Warnings, alarms and misapplications

Warning/alarm	Report	LCD display
Displayed in LCD: Battery status low	Replace battery	
Recognition function with a long beeping sound. Intensity is reset to 0 mA	Apply the adhesive electrodes to the skin	

Application programs

TENS application programs

Fixed programmed TENS application programs P01–P27

Program	Length of use	Frequency Hz	Pulse length μs
P01	30 min.	80 Hz	120 μs
P02	30 min.	100 Hz	80–140 μs
P03	30 min.	80 Hz	180 μs
P04	30 min.	100 Hz	250 μs
P05	5 min.	5 Hz	360 μs
	15 min.	1 Hz	360 μs
		40 Hz	360 μs
P06	20 min.	80 Hz	200 μs
	20 min.	100 Hz	250 μs
P07	15 min.	120 Hz	70 μs
	10 min.	100 Hz	250 μs
P08	30 min.	100 Hz	200 μs
P09	30 min.	80 Hz	150 μs
P10	30 min.	2 Hz	250 μs
P11	30 min.	100/2 Hz	150/200 μs
P12	30 min.	100/2 Hz	150/200 μs
P13	30 min.	100 Hz	150 μs
P14	30 min.	2–80 Hz	200–100 μs

P15	30 min.	80 Hz	100–180 μ s
P16	30 min.	5–15 Hz	200 μ s
P17	30 min.	80 Hz	150–200 μ s
P18	30 min.	5–15 Hz	300 μ s
P19	30 min.	80 Hz	180 μ s
P20	30 min.	10 Hz	180 μ s
P21	30 min.	100/2 Hz	150 μ s/200 μ s
P22	30 min.	2 Hz	180 μ s
P23	30 min.	80 Hz	200 μ s
P24	30 min.	2–125 Hz	100–200 μ s
P25	30 min.	80 Hz	300 μ s/180 μ s
P26	30 min.	2–120 Hz	200 μ s
P27	30 min.	100 Hz/2 Hz	150 μ s/200 μ s

Adjustable TENS programs P28, P29, P30, Pb1, Pb2 and Pb3

The values in parentheses specified in the table are the default settings. You can adjust the values in the settings (for example 5–100 Hz) in the following areas.

Pro-gram	Length of use	Frequency Hz Default setting	Pulse length μs Default setting
P28	5–100 (20) min.	1–150 (100) Hz	50–250 (200) μ s
P29	5–100 (20) min.	100 Hz	50–250 (200) μ s

P30	5–100 (20) min.	1–150 (100) Hz	100–250 μ s (pulse duration during treatment - not adjustable)
------------	-----------------	----------------	---

Program	Body area	Length of use	Frequency Hz Default setting	Pulse length μs Default setting
b1	Neck	30 min.	80 Hz	150 μ s
b2		30 min.	100 Hz	150 μ s
b3		5–100 (30) min.	20–125 (125) Hz	100–450 (200) μ s
b1	Lower back	30 min.	80 Hz/2 Hz	180 μ s
b2		30 min.	80 Hz	300–180 μ s
b3		5–100 (30) min.	20–125 (100) Hz	100–450 (330) μ s
b1	Knee	30 min.	100 Hz	200 μ s
b2		30 min.	100 Hz	150 μ s
b3		5–100 (30) min.	20–125 (80) Hz	70–150 (70) μ s

EMS application programs

Fixed EMS programs P01 to P07

Program	Length of use	Frequency Hz	Pulse length μs
P01	20 min.	50 Hz	80–300 μs
	10 min.	20 Hz	300 μs
		3 Hz	300 μs
P02	20 min.	50 Hz	60–300 μs
	20 min.	30 Hz	280 μs
		3 Hz	280 μs
P03	20 min.	30 Hz	80–260 μs
	10 min.	8 Hz	250 μs
P04	20 min.	50 Hz	60–300 μs
	20 min.	3 Hz	220 μs
		30 Hz	220 μs
P05	15 min.	55 Hz	250 μs
		4 Hz	250 μs
P06	20 min.	30 Hz	200–320 μs
	10 min.	8 Hz	400 μs
P07	20 min.	50 Hz	220–380 μs
	10 min.	3 Hz	400 μs
		20 Hz	400 μs

Adjustable TENS programs P08, P09, P10, Pb1, Pb2 and Pb3

The values in parentheses specified in the table are the default settings. You can adjust the values in the settings (for example 5–100 Hz) in the following areas.

Program	Length of use	Frequency Hz Default setting	Pulse length μs Default setting
P08	5–100 (30) min.	1–100 (30) Hz	50–320 (250) μs
P09	5–100 (30) min.	1–100 (30) Hz	100–320 μs (pulse duration during treatment - not adjustable)
P10	5–100 (30) min.	1–100 (30) Hz	(300) μs Contraction in seconds 2–30

Important:

The program P10 has a fixed pulse length of 300 μs . In programming mode, you can adjust the contraction time from 2 to 30 seconds for this program.

Program	Body area	Length of use	Frequency Hz Default setting	Pulse length μs Default setting
b1	Neck	27 min.	3–40 Hz	300 μs
b2		25 min.	3–65 Hz	300 μs
b3		5–100 (30) min.	20–125 (50) Hz	100–450 (150) μs

b1	Lower back	27 min.	3–40 Hz	250 μ s
b2		25 min.	3–65 Hz	250 μ s
b3		5–100 (30) min.	20–125 (80) Hz	100–450 (150) μ s
b1	Knee	27 min.	3–40 Hz	250 μ s
b2		25 min.	3–65 Hz	250 μ s
b3		5–100 (30) min.	20–125 (50) Hz	100–450 (350) μ s

Massage program

Fixed massage program P01 – P10

Program	Length of use	Frequency Hz	Pulse length μ s
P01	20 min.	3 Hz	320 μ s
P02	25 min.	34 Hz	80–300 μ s
P03	20 min.	25 Hz	150–300 μ s
P04	15 min.	120 Hz	80 μ s
P05	15 min.	25 Hz	80 μ s
		3 Hz	80 μ s
P06	20 min.	8 Hz	310 μ s
P07	20 min.	10 Hz	250 μ s
P08	10 min.	5 Hz	300 μ s
P09	20 min.	40 Hz	80–250 μ s
P10	20 min.	34 Hz	120–300 μ s

(Min. = minutes, Hz = cycles per second (Hertz), μ s = pulse length in microseconds)

The preinstalled programs each have a duration of max. 40 minutes. The corresponding waveforms, frequencies and pulse times for the programs are found in the program overview shown above.


Adhesive electrodes

Instructions for applying the adhesive electrodes

The following examples illustrate standard positions for the application of the adhesive electrodes. As every user responds differently to electrical nerve stimulation, for a successful pain treatment or stimulation to the muscle groups, it is very important that the adhesive electrodes are placed correctly in consultation with your physician.

Direction of electrical current

Signal form:
The signal form is bipolar pulse

- Apply the adhesive electrodes **4** above and below (or to the left and right of) the pain area.
- Do not apply the electrodes directly on the center of pain! It is important that the stimulation current is able to flow through the pain area.


Unless instructed otherwise by a physician, we recommend an average treatment time of 30 minutes, up to three times a day.

The sensation of intensity may vary from day to day. The user can adjust the intensity through the CH1/CH2 intensity buttons  .


Placement of electrodes

The illustrated examples show standard positions for the application of the adhesive electrodes. As every user responds differently to electrical nerve stimulation, for a successful pain treatment or stimulation to the muscle groups, it is very important that the adhesive electrodes are placed correctly in consultation with your physician.

Example for application of adhesive electrodes 1,57 × 1,57 in (40 × 40 mm)









Application example 21: Do not use during pregnancy.




Placement of electrodes for EMS (electrical muscle stimulation)



The placement depends on which muscle groups are to be stimulated. Various suggested positions are found in the form of pictograms as well as drawings showing the placement of the electrodes. The following chart shows different muscle groups and some useful tips on placement for the best stimulation, and how a contraction (muscle tension) can be induced.

Muscle group	Placement of electrodes	Stimulation positions	Induce contraction (muscle tension)
Sole of the foot		Sit with your feet on the floor.	Tense the muscles in your soles by trying to dig your toes into the floor.
Peroneus muscles		Sit with your feet on the floor.	Tense your calf muscles by pressing your big toe forcefully against the floor, while lifting the outside toes up from the floor at the same time.
Tibialis anterior		Sit with your feet placed under a piece of furniture, so that you can no longer bend your ankles.	Tense the shin muscles by pressing your toes forcefully upwards against an object of resistance.
Calf muscles		Sit with your back and feet supported. A simple way to do this is to sit in a doorway.	Tense the calf muscles by pressing your toes forcefully against an object of resistance.

Hamstring muscles		Lie flat on your stomach with your ankles flexed, without causing discomfort.	Tense the upper thigh muscles by attempting to raise your knee.
Adductors (inner thigh muscles)		Sit with a hard object placed between your knees (without causing discomfort).	Pull your legs to each other by trying to press your knees together.
Lower back		<p>Sitting position, please note: Due to the anatomical particularities of the lower back muscles, training in this mode requires particularly strong muscles.</p> <p>Place the electrodes on the lower back muscles as illustrated.</p>	Tense the lower back muscles by trying to sit as upright as possible.
Back muscles		Sitting position	Tense the back muscles by trying to sit as upright as possible.

Adhesive electrodes

<p>Cervical (neck) vertebrae musculature</p>	 <p>P.13</p>	<p>Sitting position</p>	<p>Tense the back muscles by trying to sit as upright as possible.</p>
<p>Trapezius muscle</p>	 <p>P.14</p>	<p>Sitting position</p>	<p>Tense the trapezius muscle by lifting and lowering the shoulders.</p>
<p>Quadriceps</p>	 <p>P.07</p> <p>Sitting position. There are two variations for this exercise: static, by blocking knee movement; or dynamic, movement against an object of resistance with the use of heavy weights.</p>	<p>Sitting position. There are two variations for this exercise: static, by blocking knee movement; or dynamic, movement against an object of resistance with the use of heavy weights.</p>	<p>Tense the quadriceps by attempting to extend your leg.</p>
<p>Buttocks</p>	 <p>P.08</p>	<p>Stand, or lie on your stomach.</p>	<p>Tense the gluteal muscles by squeezing your buttocks together.</p>

Abdominal muscles	  <p>Lie on your back, slightly elevated if desired. There are two variations for this exercise. Static: simply contracting the muscle with the movement described on the right. Dynamic: bring your torso to your thighs; in this case, take care not to place stress on the lower spine; keep the knees pressed together.</p>	<p>Tense the abdominal muscles by attempting to lift head and shoulders up from the floor.</p> <p>Attention: Do not use during pregnancy.</p>
-------------------	---	--

Technical malfunctions

Technical malfunctions and troubleshooting

Fault	Cause	Solution
The batteries have been inserted, but the display is blank.	There may be debris in the battery compartment. Check to see if the batteries are full and are inserted with the correct polarity. Check to see if the batteries are making contact.	Any foreign materials should be removed. Replace the batteries with fully charged batteries. Make sure the polarity is correct.
	There is a fault in the electronics.	Remove the batteries and put them back in after approx. 3 seconds.
The battery display is working, but the adhesive electrodes are not transmitting electric current.	The lead wires are not properly assembled.	Check that the connections to the device and to the adhesive electrodes are good.

The device is set to an intensity level, but only mild stimulation is felt through the adhesive electrodes.	The batteries are weak.	Replace the batteries with fully charged batteries. Make sure the polarity is correct.
	The surface of the skin is not clean.	Clean the skin.
	The entire adhesive surface of the electrode has no adhesive effect and is worn out.	The adhesive electrodes must be replaced.
The stimulation current becomes stronger, although the device is set to a low intensity.	The adhesive electrodes are not fully attached to the surface of the skin.	Press the adhesive electrodes firmly to the skin's surface.
	The adhesive electrodes only partially adhere to the skin.	The adhesive electrodes are old and must be replaced.
The device stops while in use.	The batteries are drained.	Replace the batteries with fully charged batteries. Make sure the polarity is correct.
	There is a fault in the electronics.	Remove the batteries and put them back in after approx. 3 seconds.
There are changes to or reddening of the skin surface.	This is possibly caused by the adhesive electrodes.	Stop use immediately and consult your physician.

Electrical interference immunity

Tips for electromagnetic interference immunity

Guidelines and manufacturer's declaration – electromagnetic emissions

This device is intended for use in electromagnetic environments as described below. The customer or operator is responsible for observing that the device is used in such an environment.

Emission test	Conformity	Electromagnetic environment – guidelines
HF emissions in accordance with CISPR 11	Group 1	The device uses RF (radio frequency) energy for its internal functions. Its HF emissions are minimal, and disturbances to nearby electronic devices are unlikely.
HF emissions in accordance with CISPR 11	Class B	Conformant
Harmonic current emissions in accordance with IEC 61000-3-2	Not applicable	This medical electrical device runs exclusively on battery power.
Voltage fluctuations/flicker emissions in accordance with IEC 61000-3-3	Not applicable	

Guidelines and manufacturer's declaration – electromagnetic interference immunity

This device is intended for use in electromagnetic environments as described below. The customer or operator is responsible for observing that the device is used in such an environment.

Interference immunity tests	Test level in accordance with IEC 60601	Degree of conformity	Electromagnetic environment
Electrostatic discharge (ESD) in accordance with IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Conformant	Floor surfaces should be wood, concrete or ceramic tile. If the floor has a synthetic covering, the relative humidity must be at least 30 %.

Electrical interference immunity

Fixed RF transmission IEC 61000-4-3	10 V/m 80 MHz 2.7 GHz	Conformant	The field intensities of fixed RF transmission equipment outside protected locations, as determined in
Conducted RF IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM and amateur radio bandwidths between 0.15 MHz and 80 MHz	Not applicable	<p>an electromagnetic location test, should not exceed 3 V/m. Interference may occur in proximity to equipment labeled as such:</p> <p> See below for calculation of the separation distance. If the transmission position is known, the specific distance can be calculated by means of available equations.</p>

Fast transient electrical bursts in accordance with IEC 61000-4-4	± 2 kV for power supply lines ±1 kV for input/output cables	Not applicable	This medical electrical device runs exclusively on battery power.
Electricity surges IEC 61000-4-5	± 0.5 kV, ± 1 kV cable(s) to cables ± 0.5 kV, ± 1 kV, ± 2 kV cable(s) to ground	Not applicable	
Power supply frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m	Conformant	Current frequency magnetic fields should correspond to the characteristics of a typical location in a typical commercial or hospital environment.

Electrical interference immunity

<p>Voltage drop, brief interruptions and voltage fluctuations in power supply lines IEC 61000-4-11</p>	<p>0 % UT; 0.5 cycle at 0°, 45 °, 90 °, 135 °, 180 °, 225 °, 270 ° and 315 °</p> <p>0 % UT; 1 cycle and 70 % UT; 25/30 cycles, single phase: at 0°</p> <p>0 % UT; 250/ 300 cycles</p>	<p>Not applicable</p>	<p>This medical electrical device runs exclusively on battery power.</p>
--	---	-----------------------	--

Calculation of separation distance for non-life-supporting systems (3Vrms / 3 V/m conformity)			
Maximum electrical output of transmission device (W)	Distance depending on frequency of transmission device (m)		
	150 kHz to 80 Mhz d=1.2	80 Mhz to 800 Mhz d=1.2	800 Mhz to 2.7 Ghz d=2.3
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.77
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

Important information on electromagnetic compatibility (EMC)

A high number of electronic devices such as PCs, mobile telephones, etc., may cause electromagnetic disturbances in the use of medical electronic devices. Electromagnetic disturbances may lead to failure in medical devices and possibly result in dangerous situations. Medical products must also not lead to faults in other devices.

The standard IEC60601-1-2 was introduced in order to meet EMC (electromagnetic compatibility) requirements and prevent dangerous situations with such devices. This standard specifies the immunity of devices to electromagnetic interference and the maximum values for electromagnetic emissions for medical devices.

Medical devices must comply with IEC60601-1-2 regarding interference immunity and emissions.

Special precautionary measures must be observed, however:

- The use of parts and wires other than those specified by the manufacturer may cause excessive emissions or reduce the interference immunity of the equipment. This excludes lead wires and replacement parts for internal components which are sold by the manufacturer.
- Medical devices should not be operated on, under, or near other devices. If operation of the device on/under/near other devices is unavoidable, the medical device must be tested for normal function in the configuration to be used.
- Information on the EMC environment in which this product should be used can be found in the EMC electric conductance table.

Cleaning and care of the TENS/EMS device

The TENS/EMS device should not be placed in direct sunlight or on hot surfaces.

Do not clean the TENS/EMS device when it is switched on and/or when the adhesive electrodes are connected.

Clean the surfaces of the TENS/EMS device gently with a soft, slightly dampened cloth. Make sure that no liquids penetrate the housing. A mild cleaning agent may be used for stubborn dirt. Clean the TENS/EMS only when it is switched off. Remove the batteries before cleaning the device. Allow the TENS/EMS device to dry thoroughly. Do not use chemical or abrasive cleaning agents for cleaning the TENS/EMS device or the adhesive electrodes.

For hygienic reasons, do not share the adhesive electrodes with other users.

Disinfect the device with a suitable commercial disinfecting agent. Allow the TENS/EMS device to dry thoroughly. Do not submerge the TENS/EMS device in water or other liquids.

Notes on using the TENS/EMS device

After each use please reattach the adhesive electrodes to protective film and put into the polybag. Please store device, adhesive electrodes and all accessories in a dry and cool place, preferably in the storage bag.

The TENS/EMS device does not require regular maintenance. Do not disassemble or attempt to repair the TENS/EMS device, as this could lead to technical mishaps or bodily injury. **Warning! Risk to life!**

Remove the batteries from the device when it will not be used for longer periods.

Information can be obtained by contacting our Service Centre.

Replacement electrodes can be ordered from our after sales service partner:

Greater Concepts By Design Team

1803 Springside Drive

Plainfield IL 60586

1-815-342-5850

greaterconcepts@live.com

Technical data, symbols, pictograms

Model:	TEN 603
Article number:	96373
Dimensions (L × W × H):	approx. 5.2 × 2.5 × 0.9 in (133 × 63 × 23 mm)
Weight:	approx. 0.16 lb (95 g) (without batteries)
Adhesive electrode area:	1.57 × 1.57 in or 2.46 in ² (40 × 40 mm or 16 cm ²)
Material:	Plastics, metals

LOT Batch name: V3318TEN603

SN Serial number: 00001
(sequential number)

 2018-04 Date of manufacture 2018-04 (year, month)

 Handelshaus Dittmann GmbH
Kissingener Strasse 68, D-97727
Fuchsstadt/Germany.

 Type BF (body float) protection against electrical shock. A type BF application equipment with greater protection against electric shock to the body, but not directly to the heart.

Device nameplate:

Electrical data

Power supply:	4.5 V DC, 3 x AAA batteries
Pulse voltage (V):	0–50 volts with a charge of 500 ohms per channel
Frequency (Hz):	TENS 1–150 Hz, EMS 1–125 Hz, MESSAGE 3–120 Hz (oscillations per second)
Pulse width (duration):	TENS 50–450 μ s, EMS 50–450 μ s, MESSAGE 80–320 μ s (microseconds)
Pulse intensity (mA):	0–90 mA with 500 ohms charge, direct current (mA= milliampere, ohm= electrical resistance)
Power consumption:	< 300 mA
Electrical tolerance:	+/- 20 % with 500 ohms charge
Output signal form:	Bipolar (symmetric-biphasic)
Application data:	2 channels with separately adjustable intensity
Ambient temperature:	Max. 41 °F – 104 °F (5 °C – 40 °C)
Air humidity with normal operation:	30 % – 75 % (percent)
Atmospheric pressure:	700 hPa – 1060 hPa (hectopascal)

Storage/transport data

 Storage/transport temperature:	-50 °F – 131 °F (-10 °C – 55 °C)
 Max. humidity for storage and transport:	10 % – 90 % (percent)
 Atmospheric pressure:	700 hPa – 1060 hPa (hectopascal)

Disposal

Disposing of the packaging

Sort the packaging before you dispose of it. Dispose of paperboard and cardboard with the waste paper and wrappers in the recyclable waste.

Disposing of old appliances

- Should the TENS/EMS device no longer be capable of being used at some point in time, dispose of it in accordance with the regulations in force in your city or state.
- Please ensure your recycling information applies to the local regulations and EPA recommendations (www.epa.gov).

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Handelshaus Dittmann GmbH Kissinger Straße 68, D-97727 Fuchsstadt/Alemania.